# GYNEX®

### Winter 2024

# **SEE** the **Difference** with

pages 46-51

# Steve says . .



As we step into the promising year of 2024, the team at Gynex extends heartfelt wishes for a Happy New Year filled with joy, well-being, and prosperity.

At Gynex, our passion for women's healthcare is the driving force behind everything we do. As we embrace the opportunities in the year ahead, we are delighted to present our latest winter catalog featuring a range of OBGYN instruments tailored to meet the evolving needs of your practice. It is our commitment to provide you with the tools that empower you to deliver the best care possible. From innovative solutions to precision instruments, each item in our catalog is a testament to

our dedication to supporting your practice. Whether you are a single office or a large hospital, we maintain a commitment to quality and value.

We recognize the importance of the work you do, and it is an honor to contribute to your success by providing reliable, top-notch instruments. Should you have any inquiries, require personalized assistance, or wish to discuss custom solutions, our team is here to ensure your satisfaction.

Thank you for choosing Gynex as your partner in women's healthcare. We look forward to serving you and being an integral part of your continued success.

Wishing you a prosperous and fulfilling year ahead!

Thank you,

Stephen J. Sullivan

Stephen J. Sullivan President and CEO, GYNEX<sup>®</sup> Corporation

- 1996 1999 International Medical Consultant
- 1993 1995 President, CooperSurgical<sup>®</sup> International Division
- 1991 1993 President, Euro-Med<sup>™</sup> Division of CooperSurgical<sup>®</sup>
- 1978 1991 Founder and President of Euro-Med<sup>®</sup>, Inc.



### **Table of Contents**

| Monsel's Ferric Subsulfate Paste                                | 2     |
|---|-------|
| Acetic Acid   | 2     |
| Reach Kit <sup>™</sup> - Cervical Block Procedure Tray          | 4     |
| Patel™, Multi-Use Forceps                                       | 5     |
| BIOPSY PUNCHES  | 6-14  |
| Ceramix <sup>™</sup> Biopsy Punch                               | 6     |
| Opal Rainbow™ Biopsy Punch, Titanium                            | 7     |
| Amethyst Biopsy Punch, Titanium                                 | 8     |
| Sapphire Biopsy Punch, Titanium                                 |       |
| Gold Biopsy Punch, Titanium                                     | 9     |
| Rose Gold Biopsy Punch, Titanium                                | 9     |
| Tischler Biopsy Punch, Stainless                                |       |
| Rotating Biopsy Punches, Stainless - ROTO Tip <sup>™</sup>      | 12-13 |
| Ring Handle Biopsy Punches, Stainless                           | 14    |
| Other Biopsy Punches, Stainless                                 | 14    |
| HOOKS   |       |
| ENDOMETRIAL CURETTES  | 15    |
| ENDOCERVICAL CURETTES   | 16    |
| Lickrish Spiral Curettes  | 17    |
| Curette Sharpening File   | 17    |
| TENACULUMS  | 18-19 |
| ENDOCERVICAL SPECULUM   |       |
| FOERSTER FORCEPS  | 19    |
| Bariatric Tissue Forceps, 1x2 teeth                             |       |
| OVUM FORCEPS  | 20-21 |
| KEYES PUNCHES   | 22    |
| Word Bartholin Catheter Kit                                     | 23    |
| Patel <sup>™</sup> , Multi-Use Forceps                          | 24    |
| IUD SECTION   |       |
| CERAMIX <sup>™</sup> Coated IUD String Scissors                 | 25    |
| IUD String Scissors   | 25    |
| Villa™ Hormonal Dissecting Forceps                              |       |
| Villa <sup>™</sup> Hormonal Implant Remover                     |       |
| IUD String Scissors, Single-Use                                 |       |
| Alligator Grasper by GYNEX®                                     |       |
| Serpentine Grasper  |       |
| IUD Removers  |       |
| Ascona® IUD Snare   |       |
| IUD Insertion Kit -GYNEX <sup>®</sup> , Single-Use <sup>®</sup> | 30    |
| UTERINE SOUNDS  |       |
| Accu-Sound®   |       |
| Sims Sound • flex   |       |
| Other Sims Sound  |       |
| UTERINE DILATORS - Stainless Steel, Soft and Dennis             |       |
| UTERINE CURETTES - Sims / Bumm                                  |       |
| Uterine Manipulator / Injector - Century UMI™                   |       |
| Endometrial Curette - SURESAMPLE® Select1                       |       |
| Endometrial Sampler - SURESAMPLE®                               |       |
| COLPOSCOPES - Galileo®  |       |
| Accessories for Galileo® Colposcope                             |       |

| Galileo® Binocular Microscope   | 53      |
|---|---------|
| COLPOSCOPES - CO Series   |         |
| Galileo® Medical Exam Light   | 60-61   |
| Accessories for Galileo® Medical Exam Light                                   |         |
| Touch-Guard® Single-Use Protective Sheath                                     |         |
| HYSTERECTOMY CLAMPS - Zephyr <sup>™</sup> Titanium Rainbow <sup>™</sup>       | 62      |
| HYSTERECTOMY SCISSORS - Titanium Rainbow <sup>™</sup> / Zephyr <sup>™</sup> . | 63      |
| NEEDLE HOLDERS - Webster / Mayo Hegar / Crile Wood                            | 64      |
| Laminaria Insertion Forceps   | 64      |
| Electrode Extender / Electrode Handpieces                                     | 65      |
| LLETZ ELECTRODES  | 66-68   |
| Electrosurgery Patient Dispersive Pads for LLETZ Procedure                    | 69      |
| Reach Kit™ - Cervical Block Procedure Tray                                    | 70      |
| Extended Needle Adapter and Accessories                                       | 71      |
| Smoke Evacuation Systems (Prefilter, Reducer, and Filter)                     | 72      |
| Tubing and Connectors for LLETZ Procedure                                     | 73      |
| Sponge Forceps, Single-Use (8)  | 74      |
| Topical Anesthetic - HurriCaine®  | 74      |
| STERILIZATION NOTICE for LLETZ Coated Instruments                             | 74      |
| LLETZ COATED INSTRUMENTS  |         |
| Tenacumlum  | 75      |
| Hooks   | 75      |
| Forceps   | 75      |
| Rulers  | 75      |
| 4- Way Expander Vaginal Speculum  | 76      |
| Lateral Vaginal Retractor   | 77      |
| Open-Sided Speculum   | 77      |
| Sullivan Vaginal Speculum   | 78-79   |
| Pederson & Graves Vaginal Speculum  | 80-81   |
| SPECULUM, STAINLESS STEEL   |         |
| Pederson & Graves Speculum  | . 82-83 |
| Sullivan Speculum   |         |
| Lateral Vaginal Retractors  | 86      |
| Auvard Weighted Speculum  | 86      |
| Guttman Speculum  | 87      |
| Greeberg Hysteroscopy Speculum  |         |
| Open-Sided Speculum   |         |
| Titanium Rainbow™ Speculum  |         |
| Fetal Doppler / Doppler Gel / Thermasonic Gel Warmer                          |         |
| pHem-Alert® Vaginal pH test   | 91      |
| Instrument Care and Cleaning  |         |
| Hot Glass Beads Sterilizer - Touch-Guard®                                     |         |
| Low Suds Detergent Concentrate - Touch-Guard®                                 |         |
| Stain Ex™ Stain Remover - Touch-Guard®  |         |
| Instrument Lubricant - Touch-Guard®   |         |
| Instrument Repair and Sharpening  |         |
|   |         |
| ORDER INFORMATION   |         |
| ORDER FORM  | 97      |

1

5% Discount on Prepaid Orders by Check or ACH payment (offer excludes microscopes, colposcopes and/or extended warranties).

# **Monsel's** Ferric Subsulfate Paste 20-22%

Monsel's is a pre-thickened formulation of the popular hemostatic agent. Apply directly to the surgical site following excision of tissue.

- Box of 12 Individual 8ml vials
- 12 rayon-tipped applicator swabs
- Single-Patient-Use (2)

**REF 5055** Box of 12 ...... **\$194**<sup>95</sup>

Helps control bleeding after biopsy, LLETZ and laser procedures.



# Acetic Acid 5%

- Box of 12 Individual 12ml vials
- 12 rayon-tipped applicator swabs
- Single-Patient-Use (2)

Helps enhance the appearance of abnormal areas.



# Galileo®

# **Family of superb optics**

Medical Exam Light see **pages 60-61** 

0 0

Colposcopes Direct-Access see **pages 46-52** 

dia.

Binocular Microscope see **page 53** 

# Reach Kit

# **For Practical Paracervical Blocks**

A convenient and cost-effective solution.

### Each Kit Contains:

- Syringe, 10ml
- Draw-Up Needle, 18 GA
- Needle Extender, 3" or 4"
- Injection Needle, 0.5" 27 GA

### Available in two size options:

| 3" Needle Exte | ender         |                 |
|----------------|---------------|-----------------|
| REF 5130       | Box of 6 Kits | \$ <b>34</b> 95 |

| 4″ | Needle  | Extender  |
|----|---------|-----------|
|    | 1100010 | EXtoniaoi |

| REF 5140 | Box of 6 Kits | \$ <b>34</b> 95 |
|----------|---------------|-----------------|
|----------|---------------|-----------------|

# **HUGE SAVINGS!**

75% less than competitor's extended needles!

To learn more, see **page 70** 



New

# **Patel**<sup>™</sup> Multi-Use Forceps



Introducing the **Patel**<sup>™</sup> Multi-Use Forceps - One instrument that every gynecologist needs in their practice.

With its patent-pending design, the **Patel**<sup>™</sup> excels at multiple functions:

- IUD Remover
- Sponge / Foerster Forceps
- Needle Holder
- Atraumatic Tissue Grasper
- More!

REF NHL3450 ...... \$114<sup>95</sup>

# The Patel is the Swiss Army Knife of Gynecology

To learn more, see **page 24** 

# **CERAMIX<sup>™</sup> Biopsy Punches**

Our sharpest, longest lasting biopsy punch


CERAMIX<sup>™</sup> coated biopsy punches offer increased resistance to corrosion, resilience against impact and low friction coefficient. Sharper up to 10 times longer than ordinary stainless steel with incomparable cutting-edge hardness. See Hardness Comparison Chart below.

- Single anchoring tooth on top and bottom jaws
- Extremely rust and corrosion resistant
- Smoother mechanical action
- Low friction coefficient
- Tip Locking Tab

### Hardness Comparison Chart

| METAL            | (HV 0.05) | COATED METAL     | (HV 0.05) |
|------------------|-----------|------------------|-----------|
| Stainless Steel  | 550       | Titanium Nitride | 2900      |
| Tungsten Carbide | 650       | CERAMIX Coating  | 4500      |

Vickers Hardness Test - HV/0.05 meaning a 50g load was used.



### **CERAMIX**<sup>™</sup>



CERAMIX<sup>™</sup> Mini-Tischler Actual Bite Size: 4.2mm x 2.3mm ■ REF 1060

20cm Shaft ......\$499<sup>95</sup>





CERAMIX<sup>™</sup> Tischler Actual Bite Size: 7mm x 3mm

**REF 1070** 20cm Shaft ......<sup>\$</sup>499<sup>95</sup>



CERAMIX<sup>™</sup> Apex<sup>™</sup>

Actual Bite Size: 7mm x 3mm

**REF 1075** 20cm Shaft ......\$499<sup>95</sup>

# **Opal Rainbow<sup>™</sup> Biopsy Punches**, TITANIUM

Designed for excellence • Manufactured for durability

Exclusive Titanium coating extends the life of instruments over traditional stainless-steel. The titanium coated biopsy punches offer increased resistance to corrosion, resilience against impact and low friction coefficient. Stays sharper and smoother longer.

- Double-spring handle
- Tip Locking Tab

GYNEX **Original!** 



### **OPAL RAINBOW<sup>™</sup> Tischler**

Actual Bite Size: 7mm x 3mm 🛛 🌒



Tischler **REF TT010** 20cm Shaft .... **\$429**<sup>95</sup>

**REF TT035** 25cm Shaft, Extra-Long ..... <sup>\$</sup>439<sup>95</sup> bariatric



Tischler, *Tip Up* **REF TT011** 20cm Shaft .... **\$429**<sup>95</sup>

REF TT036 25cm Shaft, Extra-Long ..... \$43995 bariatric



 Tischler, Tip Down

 REF TT012
 20cm Shaft ..... \$42995

REF TT037 25cm Shaft, Extra-Long ..... \$43995 bariatric

### **OPAL RAINBOW<sup>™</sup>** Mini-Tischler

Actual Bite Size: 4.2mm x 2.3mm 💼



Mini-Tischler REF TT005 20cm Shaft .... <sup>\$</sup>429<sup>95</sup>

 REF
 TT030

 25cm Shaft, Extra-Long ....
 \$43995

 bariatric



Mini-Tischler, *Tip Up* REF TT006 20cm Shaft .... <sup>\$</sup>429<sup>95</sup>

REF **TT031** 25cm Shaft, Extra-Long ..... <sup>\$</sup>439<sup>95</sup> bariatric




Mini-Tischler, *Tip Down* **REF TT007** 20cm Shaft .... <sup>\$</sup>429<sup>95</sup>

 REF
 TT032

 25cm Shaft, Extra-Long .....
 \$439<sup>95</sup>

 bariatric



# Sapphire Biopsy Punches, TITANIUM

Designed for excellence • Manufactured for durability



### Sapphire

Tischler Actual Bite Size: 7mm x 3mm REF TB010 20cm Shaft ......\$429<sup>95</sup>



### Sapphire Mini-Tischler

Actual Bite Size: 4.2mm x 2.3mm

GYNEX Original!





# Rose Gold Biopsy Punches, TITANIUM

Designed for excellence • Manufactured for durability



### Rose Gold Tischler

Ischief

Actual Bite Size: 7mm x 3mm

20cm Shaft ......\$429<sup>95</sup>



### Rose Gold

Mini-Tischler

Actual Bite Size: 4.2mm x 2.3mm

REF PINK005 20cm Shaft ......<sup>\$</sup>429<sup>95</sup> GYNEX **Original!** 

# **Tischler Biopsy Punches** STAINLESS

The Best Selling Design in The World



GYNEX Original!

| Tischler       Mini-Tischler         Image: State in the | Tischler<br>Actual Bite Size: 7mm x 3mm   | Mini-Tischler<br>Actual Bite Size: 4.2mm x 2.3mm |
|--|---|--|
| bariatric   Tischler, Tip Up   Ter 1012   Tischler, Tip Down   Tischler, Tip Down   Tischler, Tip Down   Tischler, Tip Left   Tischler, Tip Left   Tischler, Tip Left   Tischler, Tip Right  | <b>REF 1010</b> 20cm Shaft <b>\$399</b> <sup>95</sup>   |  |
| REF 1012       20cm Shaft  |   |  |
| REF 1012       20cm Shaft  | - And | No.  |
| Tischler, Tip Down       Mini-Tischler, Tip Down         RF 1014       20cm Shaft  |   |  |
| REF 1014       20cm Shaft  | <b>REF 1012</b> 20cm Shaft <b>\$399</b> <sup>95</sup>   | <b>REF 1006</b> 20cm Shaft                       |
| REF 1014       20cm Shaft  |   |  |
| Tischler, Tip Left   REF 1020   20cm Shaft   \$29995      Mini-Tischler, Tip Left Ref 1022 20cm Shaft State Tischler, Tip Right Ref 1019 20cm Shaft \$29995 Mini-Tischler, Tip Right Ref 1021 20cm Shaft \$29995   |   |  |
| Tischler, Tip Left   REF 1020   20cm Shaft   |   |  |
| REF 1020       20cm Shaft  | S   | ALE  |
| Tischler, Tip Right       Mini-Tischler, Tip Right         REF 1019       20cm Shaft   | <b>REF 1020</b> 20cm Shaft <sup>\$</sup> 299 <sup>95</sup>  |  |
| REF 1019         20cm Shaft  |   | Mini-Tischler, Tip Bight                         |
|  | REF 1019 20cm Shaft <sup>\$</sup> 299 <sup>95</sup>   |  |

### **Medium Tischler**

- Single anchoring tooth on top and bottom jaws
- Tip Locking Tab
- High Quality Stainless Steel

Available: Standard 20cm L • Extra-Long 25cm L

Actual Bite Size: 5mm x 3mm

### **Medium Tischler**

| <b>REF</b> 1008 | 20cm Shaft |  | \$ <b>39</b> | <b>9</b> 9 | Ę |
|-----------------|------------|--|--------------|------------|---|
|-----------------|------------|--|--------------|------------|---|

| REF 1033  | 25cm Shaft, Extra Long |
|-----------|------------------------|
| bariatric |                        |



### **Tischler Kevorkian**

Single anchoring tooth on top and bottom jaws

- Tip Locking Tab
- High Quality Stainless Steel

Available: Standard 20cm L • Extra-Long 25cm L

Actual Bite Size: 9mm x 3mm

### **Tischler-Kevorkian**

| <b>REF</b> 1015 | 20cm Shaft             | <b>\$399</b> 95 |
|-----------------|------------------------|-----------------|
| REF 1037        | 25cm Shaft, Extra Long | <b>\$412</b> 95 |
| bariatric       |                        |                 |



### Kevorkian


- Traditional Kevorkian bite
- Four teeth on bottom jaw
- Tip Locking Tab

### Actual Bite Size: 9mm x 3mm

 REF
 1029
 20cm Shaft
 \$399<sup>95</sup>

# **Rotating Biopsy Punches** STAINLESS

**ROTO Tip**<sup>™</sup> rotating biopsy punches offers the added convenience of changing the biopsy tip size and shape without the expense of having additional complete biopsy punches. With a single handle, you can have several tips available. Although the Standard Tip Tischler and the Mini-Tischler are the most popular, this design provides you a cost savings option to add angled tips. These tips easily rotate, allowing you to gain better positioning for a biopsy sample.



### *ROTO Tip*<sup>™</sup> Tischler

Actual Bite Size: 7mm x 3mm

### Tischler

| REF 1210  |                                 |
|-----------|---------------------------------|
| Extra Tip | <sup>\$</sup> 275 <sup>95</sup> |



| Tischler, Tip-U | lp                              |
|-----------------|---------------------------------|
| REF 1200-30     |                                 |
| Complete Punch  | <sup>\$</sup> 460 <sup>95</sup> |

| <b>REF</b> 1230 |                             |
|-----------------|-----------------------------|
| Extra Tip       | \$ <b>275</b> <sup>95</sup> |



| Tischler, Tip-Down |
|--------------------|
| REF 1200-20        |
| Complete Punch     |

| REF 1220  |                                 |
|-----------|---------------------------------|
| Extra Tip | <sup>\$</sup> 275 <sup>95</sup> |



### *ROTO Tip*<sup>™</sup> Mini-Tischler

Actual Bite Size: 4.2mm x 2.3mm

### **Mini-Tischler**

REF 1205 Extra Tip ......<sup>\$</sup>275<sup>95</sup>

| Mini-Tischler,   | Tip-Down                        |
|------------------|---------------------------------|
| REF 1200-15      |                                 |
| Complete Punch . | <sup>\$</sup> 460 <sup>95</sup> |

| <b>Mini-Tischler,</b> <i>Tip-Up</i> |                                 |
|-------------------------------------|---------------------------------|
| REF 1200-25                         |                                 |
| Complete Punch                      | <sup>\$</sup> 460 <sup>95</sup> |

| REF 1225  |                                 |
|-----------|---------------------------------|
| Extra Tip | <sup>\$</sup> 275 <sup>95</sup> |



ROTO Tip<sup>™</sup> Oval Actual Bite Size: 5mm x 3mm ●





**ROTO Tip<sup>™</sup> Kevorkian** Actual Bite Size: 9mm x 3mm

REF 1240 Extra Tip ......\$275<sup>95</sup>







### ROTO TIP<sup>™</sup> ONLY Handle

| <b>REF</b> 1200 |                                 |
|-----------------|---------------------------------|
| Handle Only     | <sup>\$</sup> 195 <sup>95</sup> |



# **Ring Handle Biopsy Punches** STAINLESS

### Sullivan Ring Handle Punch



### Kevorkian Ring Handle Punch

- Traditional Kevorkian bite
- Four teeth on bottom jaw

### Actual Bite Size: 9mm x 3mm

**REF** 1053

25cm Shaft .....\$399<sup>95</sup>



# **Other Biopsy Punches** STAINLESS



### Apex™

- V-shaped anchoring jaw on the top and bottom
- Tip Locking Tab

Available: Standard 20cm L • Extra-Long 25cm L

Actual Bite Size: 7mm x 3mm

### Apex<sup>™</sup>

| REF 1025                     | 20cm Shaft            | \$ <b>399</b> 95                |
|------------------------------|-----------------------|---------------------------------|
| REF <b>1039</b><br>bariatric | 25cm Shaft,Extra Long | <sup>\$</sup> 412 <sup>95</sup> |



# HOOKS

### **Angle Hook**

26cm Length ..... \$4995 **REF** 3225



### **2-Prong Hook**

26cm Length...... \$5995 **REF** 3227



### **Iris Hook**

**REF** 3210 26cm Length ..... <sup>\$</sup>49<sup>95</sup>



**Emmett Tenaculum** 

26cm Length......\$4995 REF 3220



# **ENDOMETRIAL CURETTES**

### **Novak Curette**

Luer lock suction attachment

Length: 25cm

| Shaft/Tip Size: 2mm REF 3302    | •<br>\$53 <sup>95</sup> |
|---------------------------------|-------------------------|
| Shaft/Tip Size: 3mm<br>REF 3303 | •<br>\$53 <sup>95</sup> |
| Shaft/Tip Size: 4mm REF 3304    | •<br>\$5395             |



### **Randall Curette**

• Luer lock suction attachment

Length: 25cm





# **ENDOCERVICAL CURETTES**

### **Kevorkian Curette** with Basket Square Tip

Basket helps hold the sample during the curettage procedure.

| Basket Size: 10.5mm x 2mm |                            |
|---------------------------|----------------------------|
| <b>REF TC3005</b>         | 0000                       |
| Titanium Rainbow          | \$ <b>79</b> <sup>95</sup> |
| <b>REF</b> 3005           |                            |
| Stainless Steel           |                            |

### Kevorkian Curette Square Tip

Basket Size: 12mm x 2mm

| REF TC3010       |                          |
|------------------|--------------------------|
| Titanium Rainbow | <b>*79</b> <sup>95</sup> |

| <b>REF 3010</b> |                                |
|-----------------|--------------------------------|
| Stainless Steel | <sup>\$</sup> 77 <sup>95</sup> |

### Dublin Curette Round Tip

This curette provides samples when aggressive scraping is necessary for hard to sample tissue.

| REF | ТСЗ | 025 |  |
|-----|-----|-----|--|
|     |     | _   |  |



**REF 3025** 

Stainless Steel .....



**Kevorkian Curette** with Basket Round Tip

Basket helps hold the sample, and the rounded tip allows easier insertion.

Basket Size: 11.5mm x 2mm

| REF TC3015       | 600                            |
|------------------|--------------------------------|
| Titanium Rainbow | <sup>\$</sup> 79 <sup>95</sup> |
| <b>REF</b> 3015  | 620                            |
| Stainless Steel  | <sup>\$</sup> 77 <sup>95</sup> |

### **Townsend Curette**


This curette has the benefit of a small pointed tip that allows easier penetration of the Os.

| REF TC3020                                  |   |                           |
|---|---|---------------------------|
| Titanium Rainbow<br>Basket Size: 12mm x 2mm |   | <b>\$79</b> <sup>95</sup> |
|   | 0 |                           |
| REF 3020<br>Stainless Steel                 |   | <b>\$77</b> 95            |
| Basket Size: 7mm x 1.3mm                    |   |                           |

### **EndoCurette ECC**

This curette is inserted into the canal and rotated to samples from the on two sides.





# **Lickrish Spiral Curettes**

The innovative Endocervical Lickrish Spiral Dual-Purpose Curette was pioneered by Dr. Gordon M Lickrish, former director of the colposcopy unit at University Health Network in Toronto. With its titanium-coated spiral design, this instrument offers enhanced resistance to corrosion and remarkable durability against impact. Specifically designed for cases of complete or severe partial cervical stenosis, this unique curette proves to be the ideal solution. It enables tissue sampling in situations where a standard curette is too large for insertion. The tip can be gently guided through the Os, while simultaneously advancing and rotating.

- Enables curettage in narrow endocervical canals, surpassing standard curette limitations
- Enables sampling in situations where a normal curette cannot access due to cervical stenosis
- Acts as a cervical dilator, ensuring precise sampling of the cervix or endometrium
- Provides dual functionality, enhancing the effectiveness and versatility of the instrument

| Titanium Cobalt<br>REF TB3040 21cm Length   | GYNCX  |
|---|--------|
| Titanium Rainbow<br>REF TC3040 21cm Length \$79 <sup>95</sup>   | GYNCX  |
| Stainless Steel<br>REF 3040 21cm Length <sup>\$</sup> 79 <sup>95</sup>  | GYNCX  |
| Curette Sharpening File<br>A small tapered diamond-coated sharpening<br>stone, helps keep your curette sharp. | GYNEX' |

### Kevorkian Curette, Single-Use

Features a consistently sharp tip modeled after the highly popular Curette with Basket and Round Tip. Designed with a firm handle that maintains stability under pressure, enabling improved maneuverability during the procedure. The angled handle provides enhanced leverage for more effective curettage.

**REF 3030** Box of 25.....**\$74**95



Total Length: 31cm Inner Basket OD: 3mm Working Length: 16cm OD: 4mm

[REF] 3030

# **TENACULUMS**

### **Braun Tenaculum**

### LiteTouch<sup>™</sup> TIP

LiteTouch tips overlap 1-2 mm ±1mm (when closed)



| Straight Handle  |             |                 |
|------------------|-------------|-----------------|
| <b>REF T2230</b> | 25cm Length | \$ <b>52</b> 95 |
| Angled Han       | dle         |                 |
| <b>REF T2235</b> | 25cm Length | \$ <b>62</b> 95 |

### **STANDARD TIP**

Standard tips overlap 4-6 mm ±1mm (when closed)



### **Straight Handle**

| REF 3502        | 25cm Length | \$ <b>52</b> 95 |
|-----------------|-------------|-----------------|
| <b>REF B230</b> | 30cm Length | \$ <b>59</b> 95 |

### Angled Handle

**REF 3532** 25cm Length ...... **\$62**<sup>95</sup>



### Schroeder Tenaculum

### LiteTouch<sup>™</sup> TIP

LiteTouch tips overlap 1-2 mm ±1mm (when closed)



| Straight Handle  |             |                 |
|------------------|-------------|-----------------|
| <b>REF T2222</b> | 25cm Length | \$ <b>52</b> 95 |
| Angled Han       | dle         |                 |
| <b>REF T2225</b> | 25cm Length | \$ <b>62</b> 95 |

### **STANDARD TIP**

Standard tips overlap 4-6 mm ±1mm (when closed)



| Straight Handle |             |                 |
|-----------------|-------------|-----------------|
| <b>REF</b> 3500 | 25cm Length | \$ <b>52</b> 95 |
| <b>REF B130</b> | 30cm Length | \$ <b>59</b> 95 |
|                 |             |                 |

### Angled Handle

| REF 3530 | 25cm Length | \$ <b>62</b> 95 |
|----------|-------------|-----------------|
|----------|-------------|-----------------|

### **DOUBLE TOOTH**



### **Straight Handle**

| REF 3501 | Len |
|----------|-----|
|----------|-----|

igth 25cm ..... **\$52**95

### **Teale Tenaculums**

### Slight Curve

REF 3536

23cm Length ...... \$5295



# ENDOCERVICAL SPECULUM

Length: 27cm

### **RATCHET HANDLE**

| <b>REF</b> 2130 | 3mm Tip | <b>*189</b> 95 |
|-----------------|---------|----------------|
| <b>REF</b> 2140 | 4mm Tip | <b>*189</b> 95 |
|                 |         |                |

### NO RATCHET

| REF 2030 | 3mm Tip <sup>\$</sup> 179 <sup>95</sup> |
|----------|---|
| REF 2040 | 4mm Tip <b>\$179</b> <sup>95</sup>      |



### **Skene Tenaculums**

Slight Curve

REF 3534 24cm Length ..... \$5295



# FOERSTER FORCEPS



### **Straight Handle**

| REF 3540        | 24cm Length        |
|-----------------|--------------------|
| <b>REF B330</b> | 30cm Length \$5995 |

### Angled Handle

| REF 3542 | 24cm Length | \$ <b>62</b> 95 |
|----------|-------------|-----------------|
|----------|-------------|-----------------|

### **Bariatric Tissue Forceps**

Hold and manipulate various tissue.

| REF BTF25 | 25cm Length <b>\$57</b> 95 |
|-----------|----------------------------|
| REF BTF30 | 30cm Length <b>\$57</b> 95 |



End view - Tip Detail

# **OVUM FORCEPS**

### **Sopher Ovum Forceps**

### **SOPHER ORIGINAL**

| Length:  | 25.5cm (10″) <i>Curved Jaw</i> 11mm |   |
|----------|-------------------------------------|---|
| REF 3544 | Ratchet Handle <sup>\$</sup> 1149   | 5 |
| REF 3546 | No Ratchet                          | 5 |

### **SOPHER STANDARD**

**REF** 3556 Batchet Handle

| Length:         | 25.5cm (10") Curved Jaw 12 | 2mm                            |
|-----------------|----------------------------|--------------------------------|
| REF 3548        | Ratchet Handle             | \$ <b>114</b> 95               |
| <b>REF</b> 3550 | No Ratchet                 | <sup>5</sup> 109 <sup>95</sup> |
| Length:         | 25.5cm (10") Curved Jaw 14 | mm                             |
|                 |                            |                                |
| REF 3552        | Ratchet Handle             |                                |

### Length: 28cm (11") Curved Jaw 12mm

| <b>REF</b> 3558 | No Ratchet | <sup>\$</sup> 109 <sup>95</sup> |
|-----------------|------------|---------------------------------|

\$11495

### Length 28cm (11") Curved Jaw 14mm

| REF 3560 | Ratchet Handle | <sup>\$</sup> 119 <sup>95</sup> |
|----------|----------------|---------------------------------|
| REF 3562 | No Ratchet     | <sup>\$</sup> 114 <sup>95</sup> |

### Length: 33cm (13") Curved Jaw 12mm

| REF 3564 | Ratchet Handle | <sup>\$</sup> 124 <sup>95</sup> |
|----------|----------------|---------------------------------|
| REF 3566 | No Ratchet     | <sup>\$</sup> 119 <sup>95</sup> |

| Length: 33cm (13") <i>Curved Jaw</i> 14mm |                |                                 |
|---|----------------|---------------------------------|
| REF 3568                                  | Ratchet Handle | <sup>\$</sup> 129 <sup>95</sup> |
| REF 3570                                  | No Ratchet     | <sup>\$</sup> 124 <sup>95</sup> |



Sopher Original "lips" protect surrounding tissue



Sopher Original Jaw Open



Sopher Standard Jaw Closed



Sopher Standard Jaw Open



No Ratchet

Ratchet



Sopher Jaw Actual Size: 11mm



Sopher Jaw Actual Size: 12mm



Actual Size: 14mm

### **Bierer Ovum Forceps**

| Length:  | 28cm (11") Straight Jaw 16m                       | n                |
|----------|---|------------------|
| REF 3572 | Ratchet Handle                                    | 14 <sup>95</sup> |
| REF 3574 | No Ratchet\$1                                     | <b>)9</b> 95     |
|          |   |                  |
| Length:  | 28cm (11") Straight Jaw 19m                       | n                |
| -        | 28cm (11") Straight Jaw 19mr<br>Ratchet Handle\$1 |                  |

| Length: 33cm (13") Straight Jaw ' | 16mm            |
|-----------------------------------|-----------------|
| REF 3580 Ratchet Handle           | <b>\$119</b> 95 |
| REF 3582 No Ratchet               | <b>\$114</b> 95 |

| Length: 33cm (13") Straight Jaw 1 | 9mm              |
|-----------------------------------|------------------|
| REF 3584 Ratchet Handle           | <b>\$124</b> 95  |
| REF 3586 No Batchet               | \$ <b>119</b> 95 |



Bierer, Jaw Open

### **Jacobson Ovum Forceps**

| Length: 2 | 28cm (11″) St | traight Jaw 20mm                |
|-----------|---------------|---------------------------------|
| REF 3588  | No Ratchet    | <sup>\$</sup> 114 <sup>95</sup> |

### Length: 33cm (13") Straight Jaw 20mm **REF** 3590 No Ratchet ......<sup>\$</sup>114<sup>95</sup>

Jacobson, Jaw Open

Blumenthal

### **Blumenthal Ovum Forceps**

| Length:  | 28cm (11″) S | traight Jaw 20mm                |
|----------|--------------|---------------------------------|
| REF 3592 | No Ratchet   | <sup>\$</sup> 114 <sup>95</sup> |



Blumenthal "lips" protect surrounding tissue

# SALE

### Saenger Ovum Forceps


| Length: 2       | 26.5cm (10.5″) Straight Jaw | <sup>,</sup> 16mm         |
|-----------------|-----------------------------|---------------------------|
| REF 3596        | Ratchet Handle              | \$ <b>114</b> 95          |
| REF 3598        | No Ratchet                  | \$ <b>109</b> 95          |
| Length:         | 26.5cm (10.5″) Curved Jaw   | 16mm                      |
| REF 3600        | Ratchet Handle              | \$ <b>114</b> 95          |
| REF 3602        | No Ratchet                  | \$ <b>109</b> 95          |
| Length:         | 29cm (11.5″) Straight Jaw ′ | 16mm                      |
| REF 3604        | Ratchet Handle              | 1 <b>29</b> <sup>95</sup> |
| <b>REF</b> 3606 | No Ratchet                  | 124 <sup>95</sup>         |
| Length:         | 29cm (11.5″) Curved Jaw 1   | 6mm                       |
| REF 3608        | Ratchet Handle              | 1 <b>29</b> <sup>95</sup> |
| REF 3610        | No Ratchet                  | 5 <b>124</b> 95           |



Saenger Inner diameter teeth within oval jaws

## **Keyes Vulvar Punch**

Features a sharpened round tip, ideal for obtaining vertical core samples from the vulvar epithelium. Designed with a sturdy handle specifically for easy ejection of lodged specimens.



### Keyes Vulvar Punch, Single-Use

STERILE EO Single-Use (Single-Use)

This Single-Use Keyes Vulvar biopsy punch has a ribbed sturdy handle for maximum control and comfort. Box of 25.



# **Keyes Vulvar Punch Set**

This essential Keyes Vulvar Punch set includes a sturdy universal handle with six high quality stainless steel interchangeable screw tips, ranging from 2mm to 7mm in diameter. The set comes in a storage case with numerically designated holding pins for each tip to stabilize tip while attaching handle. **Tip Sizes:** 2mm - 7mm



# Keyes Punch Accessories Kit

The perfect accessory kit for your Keyes Punch set. This kit includes a small forceps for lifting the tissue sample and a pair of scissors for snipping it off. The diamond-coated sharpening stone will keep your Keyes punches like new.

### **Each Kit Contains:**

- Small Pick-Up Forceps
- Small Scissors
- Diamond Coated Sharpening Stone



# Word Bartholin Catheter Kit




STERILE EO Single-Use (S) Not made with natural rubber latex.

The Word Bartholin Catheter Kit contains everything you need for an effective office treatment of Bartholin gland abscesses and cysts providing patients with immediate relief.

### Each Kit Contains:

- (1) Silicone inflatable catheter
- (1) Scalpel
- (1) 3cc Syringe w/24 GA Needle

Outpatient Treatment for Bartholin gland cysts and abscesses.





Patel<sup>™</sup> Multi-Use Forceps - One instrument that every gynecologist needs in their practice. With its patent-pending design, the Patel<sup>™</sup> excels at multiple functions!





Atraumatic Tissue Grasper



IUD Remover



Sponge Forceps



New

Needle Holder



## IUD SECTION page 25 - page 30

### **ONYX IUD String Scissors**

Blunt Tip, Slight Curve, Fine Serration

Onyx Titanium Coated Stainless Steel

### **KEY ADVANTAGES**

- Incomparable cutting-edge hardness
- Extremely rust and corrosion resistant
- Smoother mechanical action
- Low friction coefficient

Length: 25cm

REF 29-040 ..... \$104<sup>95</sup>



# *Our popular IUD String Scissors are now more durable and sharper than ever!*



### **IUD String Scissors**

Blunt Tip, Slight Curve, Fine Serration

High Quality Stainless Steel


Length: 25cm

REF 29-030 ..... \$79<sup>95</sup>



### VILLA™ New **Hormonal Implant Dissecting Forceps**


The Villa<sup>™</sup> Hormonal Implant Dissecting Forceps and Implant Remover make the perfect pair in the event that an implant is not visible upon incision. After making an incision, some implants remain difficult or impossible to view. If this is the case, the dissecting forceps, with pointed tip, can be opened and closed to separate the tissue encapsulating the implant. This will allow the remover to grasp, hold, and retract the implant.




# VILLA™ Hormonal Implant Remover





The Villa<sup>™</sup> Hormonal Implant Remover is a clamp specially designed to grasp, hold, and retract a contraceptive implant located in the arm, such as Nexplanon<sup>®</sup>. A small, ring tip secures around the mid-point of the implant for modified "U" technique.




Designed to grasp, hold, and retract a contraceptive implant located in the arm, such as Nexplanon®, using a modified "U" technique.












### **ASCONA<sup>™</sup> IUD String Snare**

### Allows gentle removal of IUD

The rotating handle provides a way to the string in the hook tip. Pulling back on the handle secures the string at the tip. Once secured, the Ascona<sup>™</sup> allows the user to gently pull out and remove the IUD. The Ascona<sup>™</sup> is easily disassembled for cleaning and sterilization. The locking tab can be loosened to remove the hook and flexible shaft.

REF 29-350 \$79<sup>95</sup>



### Gynex<sup>®</sup> IUD Insertion Kit, Single-Use

### STERILE EO Single-Use (8)

Not made with natural rubber latex.

- Reduces the risk of transmission of infectious agents
- Eliminates costly instrument sterilization

**REF 29-500** Box of 12......

### EACH KIT TRAY CONTAINS:

- (1) Tenaculum Forcep
- (1) IUD String Scissor
- (2) Rayon Swabs
- (1) ACCU-SOUND® Uterine Sound
- (1) Folded Drape




# **UTERINE SOUNDS**

### ACCU-SOUND® Uterine Sound, Single-Use

STERILE EO Single-Use S

The simple yet ingenious design of the ACCU-SOUND<sup>®</sup> Uterine Sound eliminates the guesswork involved with a standard uterine sound. When inserting the ACCU-SOUND<sup>®</sup> into the uterus, the ring rests against the outside of the Os. When the device is withdrawn from the uterus and vagina, the ring stays in place, clearly indicating the uterine measurement.



| Sims Sound • flex, Single-Use   | New  |
|---|--|
| STERILE EO Single-Use ®   | GYNEX<br>Original!   |
| <ul> <li>Gynex has improved the Sims Sound. Engineered from the Sims Sound Flex formable shaft enhances the proveach patient's unique physical requirements. Bold, lase clear, reliable measurements.</li> <li>Formable Memory Polymer</li> </ul> | iders ability to pre-shape the sound to meet                             |
| • OD 4mm  |  |
| Working Length: 16.6cm<br>REF 4630 Box of 25\$54 <sup>95</sup>  |  |
| a to   | I* I* GYNEX*   |
|   | Not Shown at Actual Size   |
| Other Sims Sounds<br>Sims Sound Plus+, Single-Use   | STERILE EO Single-Use (S)<br>Not made with natural rubber latex.         |
|   |  |
| Malleable plastic     OD 4mm  | BEFERERE GYNEX   |
| • OD 4mm  |  |
|   | GYNEX  |
| • OD 4mm<br>Length: 24cm  | GYNEX  |
| • OD 4mm<br>Length: 24cm<br>REF 4675 Box of 25\$5695  | STERILE EO Single-Use S<br>Not made with natural rubber latex.           |
| <ul> <li>OD 4mm</li> <li>Length: 24cm</li> <li>REF 4675 Box of 25<sup>\$5695</sup></li> <li>Sims Sound, Single-Use</li> <li>Malleable plastic</li> <li>OD 4mm</li> <li>Length: 16.6cm</li> </ul>  | STERLE EO       Single-Use ®         Not made with natural rubber latex. |
| <ul> <li>OD 4mm</li> <li>Length: 24cm</li> <li>REF 4675 Box of 25</li></ul>   | STERILE EO Single-Use S<br>Not made with natural rubber latex.           |

Length: 32cm

| <b>REF</b> 4334 | Stainless Steel, Non-Malleable   | \$ <b>42</b> 95 |
|-----------------|----------------------------------|-----------------|
| REF 4336        | Stainless Steel, Malleable ····· | SA <b>7</b> 95  |



# **UTERINE DILATORS**

### Mini Dilator Set

### Set of 5 inter-changeable Probes and 1 Handle

High Quality Stainless Steel **Sizes:** 1mm, 1.5mm, 2mm, 2.5mm and 3mm





| Extra | Mini | Dilator | Tips |  |
|-------|------|---------|------|--|
|-------|------|---------|------|--|

(Fits REF 4005 Handle)

**Sizes:** 4mm and 5mm Available only from GYNEX


| REF 4005.4 | 4mm | \$ <b>42</b> 95 |
|------------|-----|-----------------|
| REF 4005.5 | 5mm | \$ <b>42</b> 95 |



Dots show actual diameter in mm







### Soft Dilators, Single-Use

STERILE EO Single-Use ® Not made with natural rubber latex.

Designed to dilate the cervical canal to allow entry into the uterus for evaluation of the cervical canal or uterine cavity. These color coded dilators are soft and malleable allowing the operating physician to dilate with excellent control.

#### Soft Dilator, Os Locator

Used for locating the external Os by gently pushing the tip against the cervix.

Tapered Tip Length: 23mm

Overall Length: 21.5cm

| REF 4211 | Box of 12 | \$ <b>52</b> 95 |
|----------|-----------|-----------------|
| REF 4221 | Box of 24 | <b>\$94</b> 95  |

#### Soft Dilator, Canal Finder

Used for indicating the cervical canal.

Tapered Tip Length: 45mm

Overall Length: 21.5cm

| REF 4212 | Box of 12 | \$ <b>52</b> 95 |
|----------|-----------|-----------------|
| REF 4222 | Box of 24 | <b>\$94</b> 95  |

#### **Soft Dilator, Fundus Finder**

Used for establishing the entire length of the cervical canal from external Os to the internal Os.

Tapered Tip Length: 85mm

Overall Length: 21.5cm

| REF 4213 | Box of 12 | \$ <b>52</b> 95 |
|----------|-----------|-----------------|
| REF 4223 | Box of 24 | <b>\$94</b> 95  |





36

www.gynex.com • 1.888.486.4644



# Soft Teflon Dilators, REUSABLE

Designed to dilate the cervical canal to allow entry into the uterus for evaluation of the cervical canal or uterine cavity. These Soft Teflon Dilators are reusable, fully autoclavable and can be sterilized after use. Each of these three dilators are comfortable to hold and allow the operating physician to dilate with excellent control.

# $\checkmark$

**STERILE / Reusable** 

**OK to Autoclave** 

#### Soft Teflon Dilator, Os Locator

Tapered Tip Length: 23mm Overall Length: 21.5cm

REF 4201 \$4495

#### Soft Teflon Dilator, Canal Finder

 Tapered Tip Length: 45mm

 Overall Length: 21.5cm

 REF
 4202

#### Soft Teflon Dilator, Fundus Finder

 Tapered Tip Length: 45mm

 Overall Length: 21.5cm

 REF
 4203

# **Great Value!** SAVE on a SET of 3 Dilators

STERILE / Reusable OK to Autoclave

# Soft Teflon Dilators Set, REUSABLE

**Set of three includes:** one Os Locator, one Canal Finder and one Fundus Finder.

REF 4200

Set of 3 ..... \$9995

### Soft Teflon Dilators Set, REUSABLE NON-STERILE

**Set of three includes:** one Os Locator, one Canal Finder and one Fundus Finder.

REF 4300

Set of 3 ..... **\$94**95

Note: Multi-Colored Single-Use  $\circledast$  and the White Reusable Soft Dilators are the same shape and are intended for identical use.



# Dennis Dilators® Set, REUSABLE

#### STERILE EO

Not made with natural rubber latex.

5mm to 14mm

- Dual sizes on each dilator
- Soft plastic
- STERILE Sealed in procedure blister trays
- Reusable sterilize after use (fully autoclavable)

#### Overall length: 28cm

# Great Value!

### SAVE on a SET of 5 Dilators

**REF 4230** Set of 5 ..... **\$94**95

#### **Procedure Blister Tray**

- Safe and Secure
- Sterile
- Easy Open Cover



# **UTERINE CURETTES**

### **Sims Uterine Curette**

#### Sharp & Rigid

All Sims Uterine Curette Tips Shown at Actual Size

|   |            | $\nabla$        | $\square$ | $\bigcap$                          | $\bigcap$   | $\bigcap$                      |
|---|------------|-----------------|-----------|------------------------------------|-------------|--------------------------------|
|   | 1          | Y               | Y         | Y                                  | V           | V                              |
| Size 00 Size 0  | Size 1     | Size 2          | Size 3    | Size 4                             | Size 5      | Size 6                         |
| <b>REF CSU3900-9</b><br><b>Size 00:</b> 10.5 x 4.0mm, Sha | arp, Rigid | <b>\$68</b> 95  |           | <b>sU3900-3</b><br>: 16 x 7mm, Sha | arp, Rigid  | <b>\$68</b> 95                 |
| REF CSU3900-0<br>Size 0: 10.5 x 4.5mm, Shar               |            |                 | REF C     | <b>SU3900-4</b><br>: 20 x 9mm, Sh  |             | \$0.095                        |
| REF CSU3900-1<br>Size 1: 13 x 5mm, Sharp, F               | Rigid      | \$ <b>68</b> 95 |           | <b>SU3900-5</b><br>: 22 x 9mm, Sha | arp, Rigid  | <b>\$68</b> 95                 |
| REF CSU3900-2<br>Size 2: 16 x 6mm, Sharp, F               | Rigid      | \$ <b>68</b> 95 |           | <b>5U3900-6</b><br>: 26 x 11mm, Sł | narp, Rigid | <sup>\$</sup> 68 <sup>95</sup> |

# **Bumm Uterine Curette**



40

www.gynex.com • 1.888.486.4644

# Century UMI<sup>™</sup> Uterine Manipulator/Injector



STERILE EO Single-Use (Single-Use (Single-

The Century UMI<sup>™</sup> Uterine Manipulator is easy-to-use and designed for patient safety. Optimized to support a broad procedural versatility: Laparoscopic Tubal Occlusion, Diagnostic Laparoscopy, Minilaparotomy, Salpingoplasty and Fertility Exams.

#### Specifications

- Cannula length: 33cm
- OD 4.5mm
- Luer Lock Syringe 10cc

Now including a 10cc luer lock syringe for your convenience!

#### Features

- Pear-shaped balloon for optimal hold
- Incorporated cervical stop (no spring!)
- Reduced cannula to allow easy insertion
- Balloon forms an excellent seal at internal Os

# **NEW LOW PRICE!**

Low unit price, \$18.33 each

**REF** 10-160

Box of 12 ..... \$219<sup>95</sup>

# Passing Savings To You!

Note: It is highly recommended that the practitioner utilize a sound to determine the patient's uterine depth prior to using any uterine manipulator/injector.

# suresample • Select 1

# **Endometrial Curette**

STERILE EO Single-Use (S) Not made with natural rubber latex.

Suresample<sup>®</sup> *Select1* is a uniquely designed curette for sampling the uterine endometrium, ideal for obtaining samples from the uterus and mucosal lining by scraping or vacuum suction.

#### Features

- Randall-Style sharp cutting-edge
- Semi flexible curette with slight curve
- Highly visible 1cm markings on cannula
- Cannula length: 8-1/8" (20.5cm)
- White Cannula Stylet insert for additional rigidity
- Depth Marker acts as a visual aid to protect against perforation
- OD 3mm
- 10mL standard syringe with Twist-and-Lock plunger



| <b>REF</b> 3130 | Box of 5  | \$ <b>32</b> 95 |
|-----------------|-----------|-----------------|
| REF 3135        | Box of 24 | <b>\$134</b> 75 |

*The unique design and special features allow for better accuracy and a more reliable sample.* 

# SUR Sample • Select 1 Endometrial Curette

# Simple, Safe and Effective

Tapered tip for easier insertion. Randall - Style sharp cutting edge.

#### **DEPTH MARKER**

Set the desired curette marker depth after sounding the uterus. This depth marker acts as a visual aid to reduce the risk of perforation.

#### TWIST AND LOCK PLUNGER

The twist-and-lock plunger on the syringe allows the provider to maintain a constant vacuum pressure while collecting the specimen

# The endometrial sampler with a proven track record!

# SUTE Endometrial Sampler


The Suresample<sup>®</sup> Endometrial Sampler has been trusted by doctors for over *25 years* 

Our patented 3-port sampling system doubles the harvested samples of tested competitors and allows easy retrieval with reduced cell damage. A much greater histological adequacy and integrity than any other samplers. 12 to 1 ratio of histological adequacy and 8 to 3 for histological integrity, Suresample<sup>®</sup> is the finest endometrial sampler.

Suresample<sup>®</sup> formable shaft memory allows physicians to **pre-shape shaft** to meet each patient unique physical requirements. Its tensile strength and rigidity make it malleable enough to conform to the direction of the cervix and lower segment, but strong enough to go through a tight canal. The formable 3mm shaft diameter enhances ease of use and minimizing patient discomfort.



AMOUNT OF BIOPSY TISSUE



Better performance than other brands

Suresample<sup>®</sup> Endometrial Sampler retrieves more than twice the harvest of its competition!

# <u>sure</u>sample<sup>®</sup> **Endometrial Sampler**



#### Features and Benefits

- 3-Port Sampling System one port on each side of the distal end of the shaft and one on the tip end.
- Grip retraction creates a strong internal vacuum.
- Two phase action between the ports and grip produces a higher yield of tissue sample in an amazingly controlled manner.
- Formable memory shaft allows you to shape the sampler to meet each patient's physical requirements.
- Tensile strength and rigidity allow conformity but also provide enough strength to fit through a smaller canal.
- No cutting or splitting of tip (a tedious process required by other widely used samplers that leads to damaged specimen).
- Quick retrieval process simply push the SURESAMPLE® Endometrial Sampler grip back, forcing the • specimen out of the open sampler tip and sides.

#### **Specifications**

- STERILE, Single-Use
- 3mm shaft diameter
- 24cm shaft length

Box of 25.....\$119<sup>95</sup> **REF** 3115





# GalileO<sup>®</sup> Colposcope Direct-Access<sup>™</sup> 530VCDA

#### Three-Step Magnification with Integrated LED Light and HD Video Camera

Galileo<sup>®</sup> Colposcope Direct-Access<sup>™</sup> 530VCDA features a three-step magnification with integrated LED light and HD video camera. The offset center post provides a direct visual and easy access. Combining modern technology with classic aesthetics, the Galileo<sup>®</sup> Colposcope Direct-Access<sup>™</sup> 530VCDA is a piece of art designed to facilitate comprehensive cervical examination, with ease and style.

#### Features

- Direct-Access<sup>™</sup> Offset Center Post (Eliminates an obscured view and provides a direct visual and easy access)
- Multi-Coated, Hand-Polished Apochromatic Lens (Minimizes secondary color and corrects chromatic abrasion)
- Integrated Power Source and Cables (Enclosed in base and post for a clean and clutter free look)
- Low voltage DC power (Eliminates the need to replace fuses)
- LED Illumination System (Accurate color balance, cool to the touch, long bulb life)
- Three-Step Magnification
- Optical Measurement Rings
- Individual Eye-Shields
- Integrated Green Filter
- HDMI Video
  - (Allows for live feed and image/video capture)

#### Options

- 360° Tilt Stand Base / 5-Wheel Rolling Stand Base
- Additional Extended Warranty

#### Specifications

- Focal Length: ± 300mm
- Magnification: 3.75X, 7.5X, 15X
- Field Of View: 79mm, 39mm, 19mm
- Depth Of Field: 4.5mm, 1.13mm, 0.76mm
- Light Intensity: > 30,000 LUX
- Working Height: 975mm to 1135mm
- Power Supply: 100V to 240V
- Video Feed: HDMI-Video


46

• Warranty: 1 Year - Parts and Labor (excludes shipping charges)

Galileo<sup>®</sup> Colposcope Direct-Access<sup>™</sup> 530VCDA • Three-Step Magnification with Integrated LED Light and HD Video Camera

| with a 5-wheel Rolling Base: REF 530VCDA-R                      | <sup>\$</sup> 7,395 <sup>00</sup>                   |
|---|---|
| with a 360° rotating Tilt Base: REF 530VCDA-T                   | <sup>\$</sup> 7,395 <sup>00</sup>                   |
| 1 Year Extended Warranty: REF 530VCDA T/R-1 \$375 <sup>00</sup> | 2 Year Extended Warranty: REF 530VCDA T/R-2 \$53500 |





# GalileO<sup>®</sup> Colposcope Direct-Access<sup>™</sup> 400DA

#### Three-Step Magnification with Integrated LED Light

Galileo<sup>®</sup> Colposcope Direct-Access<sup>™</sup> 400DA features a three-step magnification with integrated LED light. The offset center post provides a direct visual and easy access. Combining modern technology with classic aesthetics, the Galileo<sup>®</sup> Colposcope Direct-Access<sup>™</sup> 400DA is a piece of art designed to facilitate comprehensive cervical examination, with ease and style.

#### Features

- Direct-Access<sup>™</sup> Offset Center Post (Eliminates an obscured view and provides a direct visual and easy access)
- Multi-Coated, Hand-Polished Apochromatic Lens (Minimizes secondary color and corrects chromatic abrasion)
- Integrated Power Source and Cables (Enclosed in base and post for a clean and clutter free look)
- Low voltage DC power (Eliminates the need to replace fuses)
- LED Illumination System (Accurate color balance, cool to the touch, long bulb life)
- Three-Step Magnification
- Optical Measurement Rings
- Individual Eye-Shields
- Integrated Green Filter

#### Options

- 360° Tilt Stand Base / 5-Wheel Rolling Stand Base
- Additional Extended Warranty


#### Specifications

- Focal Length: ± 300mm
- Magnification: 3.75X, 7.5X, 15X
- Field Of View: 79mm, 39mm, 19mm
- Depth Of Field: 4.5mm, 1.13mm, 0.76mm
- Light Intensity: > 30,000 LUX
- Working Height: 975mm to 1135mm
- Power Supply: 100V to 240V
- Warranty: 1 Year Parts and Labor (excludes shipping charges)



#### Galileo<sup>®</sup> Colposcope Direct-Access<sup>™</sup> 400DA • Three-Step Magnification with Integrated LED Light

| with a 5-wheel Rolling Base: REF 400DA-R         | \$5,795 <sup>00</sup>                             |
|--|---|
| with a 360° rotating Tilt Base: REF 400DA-T      | <sup>\$</sup> 5,795 <sup>00</sup>                 |
| 1 Year Extended Warranty: REF 400DA T/R-1\$37500 | 2 Year Extended Warranty: REF 400DA T/R-2 \$53500 |



# GalileO<sup>®</sup> Colposcope Direct-Access<sup>™</sup> 200DA

#### Single Magnification with Integrated LED Light

Galileo<sup>®</sup> Colposcope Direct-Access<sup>™</sup> 200DA features single magnification with integrated LED light. The offset center post provides a direct visual and easy access. Combining modern technology with classic aesthetics, the Galileo<sup>®</sup> Colposcope Direct-Access<sup>™</sup> 200DA is a piece of art designed to facilitate comprehensive cervical examination, with ease and style.

Multiple Paten

#### Features

- Direct-Access<sup>™</sup> Offset Center Post (Eliminates an obscured view and provides a direct visual and easy access)
- Multi-Coated, Hand-Polished Apochromatic Lens (Minimizes secondary color and corrects chromatic abrasion)
- Integrated Power Source and Cables (Enclosed in base and post for a clean and clutter free look)
- Low voltage DC power (Eliminates the need to replace fuses)
- LED Illumination System (Accurate color balance, cool to the touch, long bulb life)
- Single Magnification
- Optical Measurement Rings
- Individual Eye-Shields
- Integrated Green Filter

#### Options


- 360° Tilt Stand Base / 5-Wheel Rolling Stand Base
- Additional Extended Warranty

#### Specifications

- Focal Length: ± 300mm
- Magnification: 7.5X
- Field Of View: 39mm
- Depth Of Field: 1.13mm
- Light Intensity: > 30,000 LUX
- Working Height: 975mm to 1135mm
- Power Supply: 100V to 240V
- Warranty: 1 Year Parts and Labor (excludes shipping charges)

#### Galileo<sup>®</sup> Colposcope Direct-Access<sup>™</sup> 200DA • Single Magnification with Integrated LED Light

| with a 5-wheel Rolling Base: REF 200DA-R          | <sup>\$</sup> 4,795 <sup>00</sup>                 |
|---|---|
| with a 360° rotating Tilt Base: REF 200DA-T       |   |
| 1 Year Extended Warranty: REF 200DA T/R-1 \$37500 | 2 Year Extended Warranty: REF 200DA T/R-2 \$53500 |



# **ACCESSORIES**

# for GalileO<sup>®</sup> Colposcope Direct-Access<sup>™</sup> 530VCDA

#### **Video Capture Device**

Connects Galileo<sup>®</sup> Colposcope 530VCDA to a PC for viewing and capturing images. Uses a USB 2.0 or 3.0 port on a computer.

REF VC-Magewell .....\$29900



# Cut the Cord!

Equipped with rechargeable battery.

# **Gynex®** Monitor

The Gynex<sup>®</sup> monitor is equipped with a built-in rechargeable battery for your convenience, and can easily be recharged. It offers a convenient viewing option for patients and staff with one single cable connecting colposcope to the monitor.

Monitor with Rolling Stand Base REF 450MS ......\$2,29500

Monitor Only (without base) REF 450M ...... \$29500

(Note: Actual monitor shipped may be different than model shown)

# Galileo® Binocular Microscope

#### **Specifications**

- Condenser Numerical Aperture: Na = 1.25
- Stage Movement Range: Longitudinal 50mm, Transverse 75mm
- Fine Focusing: Minimum Division .002mm
- Interpupillary Distance Adjustment Range: 53–75mm
- Light Source: 6v / 20w Halogen Lamp, Adjustable Brightness

Extended Warranties Available (Excludes shipping charges)

# Clearer images and wider views at a great price!

| Eyepieces            | 10X Magnification         | 25mm Focus | Φ 18mm Field |      |            |
|----------------------|---------------------------|------------|--------------|------|------------|
| Semi-Plan Achromatic |                           |            |              |      |            |
|                      | Magnification             | 4X         | 10X          | 40X  | 100X (oil) |
| Objectives           | Total Magnification       | 40X        | 100X         | 400X | 1000X      |
|                      | N.A. (Numerical Aperture) | .1         | .25          | .65  | 1.25       |
|                      | W.D. (mm)                 | 23.1       | 4.1          | .6   | .38        |



#### Components

- Plastic Carrying Case
- 2 Eyepieces
- 1 Extra Fuse
- 1 Extra Bulb
- 1 Dust Cover
- 1 Instructions for Use
- 1 Bottle of Immersion Oil

Warranty: 1 Year - Parts and Labor

 Galileo® Binocular Microscope
 REF
 MI6000
 \$92495

 1 Year Extended Warranty:
 REF
 MI6000 -1
 2 Year Extended Warranty:
 REF
 MI6000 - 2
 \$17500

# CO4100 Colposcope

#### Three-Step Magnification, 45° Angle with Integrated LED Light

The CO4100 Colposcope features a three-step magnification and 45° viewing angle with integrated LED light. The 45° angle lenses allow inclined viewing for a natural sitting position, reducing the need to bend over. The CO4100 Colposcope includes excellent optics designed to facilitate comprehensive cervical examination, with ease and style.

#### **Features**

- Multi-Coated, Hand-Polished Apochromatic Lens (Minimizes secondary color and corrects chromatic abrasion)
- Integrated Power Source and Cables (Enclosed in base and post for a clean and clutter free look)
- Low voltage DC power (Eliminates the need to replace fuses)
- LED Illumination System (Accurate color balance, cool to the touch, long bulb life)
- Three-Step Magnification
- Optical Measurement Rings
- Individual Eye-Shields
- Integrated Green Filter
- Adjustable Tension Ball-Joint for Easy 360° Motion

#### **Specifications**

- Focal Length: ± 316mm
- Magnification: 3.75X, 7.5X, 15X
- Optical Angle: 45°
- Field Of View: 9.5mm, 19mm, 38mm
- Depth Of Field: 5mm, 8mm, 32mm
- Light Intensity: > 25,000LUX
- Working Height: 838mm to 1194mm (33" to 47")
- Power Supply: 90V to 240V
- Warranty: 1 Year Parts and Labor (excludes shipping charges)

#### Options

Additional Extended Warranty

**CO4100 Colposcope** • Three-Step Magnification 45° Angle with Integrated LED Light

360° rotating Tilt Base: REF CO4100-LED ..... \$4,99500

1 Year Extended Warranty: REF CO4100-1 ...... \$325°0 2 Year Extended Warranty: REF CO4100-2 ...... \$475°0



# CO4100 Colposcope

Our colposcopes now utilize low-voltage DC power eliminating the need to replace fuses!





LED Light Source

Lasts Longer Cool To The Touch Color Balance



#### **Universal Teflon Ball Joint**

360° Easy Maneuver

# CO3100 Colposcope

#### Three-Step Magnification with Integrated LED Light

The CO3100 Colposcope features a three-step magnification with integrated LED light. The device is compact, and light-weight with 360° tilt base for easy maneuverability. The CO3100 Colposcope includes excellent optics designed to facilitate comprehensive cervical examination, with ease and style.

#### Features

- Multi-Coated, Hand-Polished Apochromatic Lens (Minimizes secondary color and corrects chromatic abrasion)
- Integrated Power Source and Cables (Enclosed in base and post for a clean and clutter free look)
- Low voltage DC power (Eliminates the need to replace fuses)
- LED Illumination System (Accurate color balance, cool to the touch, long bulb life)
- Three-Step Magnification
- Optical Measurement Rings
- Individual Eye-Shields
- Integrated Green Filter
- Adjustable Tension Ball-Joint for Easy 360° Motion

#### Specifications

- Focal Length: ± 316mm
- Magnification: 3.75X, 7.5X, 15X
- Field Of View: 9.5mm, 19mm, 38mm
- Depth Of Field: 5mm, 8mm, 34mm
- Light Intensity: > 25,000LUX
- Working Height: 838mm to 1194mm (33" to 47")
- Power Supply: 90V to 240V
- Warranty: 1 Year Parts and Labor (excludes shipping charges)

#### Options

Additional Extended Warranty



1 Year Extended Warranty: REF CO3100-1 ...... \$32500 2 Year Extended Warranty: REF CO3100-2 ...... \$47500



# CO3100 Colposcope

Our colposcopes now utilize low-voltage DC power eliminating the need to replace fuses!



#### LED Light Source

- Lasts Longer
- Cool To The Touch
- Color Balance



#### **Universal Teflon Ball Joint**

• 360° Easy Maneuver

# **CO1100 Colposcope**

#### Single Magnification with Integrated LED Light

The CO1100 Colposcope features 9x magnification. The device is compact and light-weight with a with 360° tilt base for easy maneuverability. The CO1100 Colposcope includes excellent optics designed to facilitate comprehensive cervical examination, with ease and style.

#### Features

- Multi-Coated, Hand-Polished Apochromatic Lens (Minimizes secondary color and corrects chromatic abrasion)
- Integrated Power Source and Cables (Enclosed in base and post for a clean and clutter free look)
- Low voltage DC power (Eliminates the need to replace fuses)
- LED Illumination System (Accurate color balance, cool to the touch, long bulb life)
- Single Magnification
- Optical Measurement Rings
- Individual Eye-Shields
- Integrated Green Filter
- Adjustable Tension Ball-Joint for Easy 360° Motion

#### Specifications

- Focal Length: ± 316mm
- Magnification: 9X
- Field Of View: 22mm
- Depth Of Field: 8mm
- Light Intensity: > 25,000LUX
- Working Height: 838mm to 1194mm (33" to 47")
- Power Supply: 90V to 240V
- Warranty: 1 Year Parts and Labor (excludes shipping charges)

#### Options

Additional Extended Warranty



58



# CO1100 Colposcope

Our colposcopes now utilize low-voltage DC power eliminating the need to replace fuses!





#### LED Light Source

- Lasts Longer
- Cool To The Touch
- Color Balance



#### Universal Teflon Ball Joint

• 360° Easy Maneuver

# **Galileo<sup>®</sup> Medical Exam Light**

Experience the ultimate convenience and functionality with the Galileo<sup>®</sup> Medical Exam Light. Its **flexible neck** and adjustable design allow you to position the light exactly where you need it, providing optimal sharp illumination. Its **durable and aesthetically pleasing** design makes it a smart investment for any medical facility.

With energy-efficient LED bulb, the Galileo® Medical Exam Light is not only cost-effective but also long-lasting.

Ideal for hospitals, clinics, and private practices, the Galileo<sup>®</sup> Medical Exam Light is the perfect addition to any medical setting.

#### **Brighter and Cost Effective**

- Features 1 white 3W LED for precise, high-output lighting
- Boasts a long lifetime of up to 50,000 hours and provides high illuminance of 70,000 Lux at a 200mm working distance
- Consumes less energy than traditional lighting options ٠
- Eliminate the need for frequent bulb replacement ٠

#### **True-to-Life Color Reproduction**

- The lighting is easy on the eyes, providing a comfortable visual experience
- Delivers exceptional color rendering capabilities, ensuring that the illuminated area looks natural and highly contrasted

#### Adjustable Spot Size

- Offers adjustable spot size with a diameter ranging from 15mm to 220mm, making it suitable for a wide range of situations and working distances, from 200mm to 1000mm
- Easy to clean, ensuring hassle-free maintenance

#### Flexible Goose Neck

- An extremely flexible goose neck can be securely fixed in any position
- Features a two-segment design that allows for bending at any angle and any directions

# **Medical Exam Light Accessories**

### **Table Mount**

Mount the Galileo<sup>®</sup> Medical Exam Light to your table for additional floor space.

REF MEL-20 ...... \$87<sup>95</sup>

#### Wall Mount

The Galileo<sup>®</sup> Medical Exam Light Wall Mount Bracket frees up floor space when necessary.

.....<sup>\$</sup>72<sup>95</sup> REF MEL-10











# Invest in the best and choose Galileo®

REF MEL-50 \$499<sup>95</sup>



SCAN to SHOP

### Touch-Guard® Single-Use **Protective Sheath**

Touch-Guard<sup>®</sup> Single-Use Sheaths protect your medical exam light and reduce the risk of contamination.

- Recommended to use with Galileo<sup>®</sup> Medical • Exam Light **REF MEL-50**
- Can be used with most goose-neck exam lights • under 4 cm in diameter





# **HYSTERECTOMY CLAMPS**



# Zephyr<sup>™</sup>Titanium Rainbow<sup>™</sup>

**Titanium Coated Hysterectomy Clamps** 

#### Straight

| <b>REF TT580-21</b> | 21cm <b>\$194</b> 95           |
|---------------------|--------------------------------|
| <b>REF TT580-24</b> | 24cm\$ <b>204<sup>95</sup></b> |
| <b>REF TT580-27</b> | 27cm\$ <b>214<sup>95</sup></b> |
| <b>REF TT580-30</b> | 30cm <b>\$224</b> 95           |
| <b>REF TT580-36</b> | 36cm <b>\$234</b> 95           |

#### **Slight Curve**

| <b>REF TT585-21</b> | 21cm <b>\$194</b> 95 |
|---------------------|----------------------|
| <b>REF TT585-24</b> | 24cm <b>\$204</b> 95 |
| <b>REF TT585-27</b> | 27cm\$ <b>214</b> 95 |
| <b>REF TT585-30</b> | 30cm <b>\$224</b> 95 |
| <b>REF TT585-36</b> | 36cm <b>\$234</b> 95 |

#### **Strong Curve**

| <b>REF TT590-21</b> | 21cm | \$ <b>194</b> <sup>95</sup> |
|---------------------|------|-----------------------------|
| <b>REF TT590-24</b> | 24cm | \$ <b>204</b> <sup>95</sup> |
| <b>REF TT590-27</b> | 27cm | \$ <b>214</b> 95            |
| <b>REF TT590-30</b> | 30cm | \$ <b>224</b> <sup>95</sup> |
| <b>REF TT590-36</b> | 36cm | \$ <b>234</b> 95            |



#### Angled

| <b>REF TT595-21</b> | 21cm | \$ <b>194</b> 95            |
|---------------------|------|-----------------------------|
| <b>REF TT595-24</b> | 24cm | \$ <b>204</b> <sup>95</sup> |
| <b>REF TT595-27</b> | 27cm | \$ <b>214</b> 95            |
| <b>REF</b> TT595-30 | 30cm | \$ <b>224</b> <sup>95</sup> |
| <b>REF TT595-36</b> | 36cm | \$ <b>234</b> <sup>95</sup> |

# The Original Zephyr<sup>™</sup>

#### Stainless Steel Hysterectomy Clamps

| Zimm               | Straigh                        | t |
|--------------------|--------------------------------|---|
| <b>REF</b> 3580-21 | 21cm <b>\$114</b> 9            | 5 |
| <b>REF</b> 3580-24 | 24cm <b>\$119<sup>9!</sup></b> | 5 |
| <b>REF</b> 3580-27 | 27cm <b>\$124</b> 9            | 5 |
| <b>REF</b> 3580-30 | 30cm <b>\$129</b> 9            | 5 |



### **Slight Curve**

| REF 3585-21 | 21cm <b>\$114</b> 95            |
|-------------|---------------------------------|
| REF 3585-24 | 24cm <b>\$119<sup>95</sup></b>  |
| REF 3585-27 | 27cm <b>\$124</b> 95            |
| REF 3585-30 | 30cm <b>\$129</b> <sup>95</sup> |



# Strong Curve

| REF 3590-21 | 21cm <b>\$114</b> 95           |
|-------------|--------------------------------|
| REF 3590-24 | 24cm <b>\$119<sup>95</sup></b> |
| REF 3590-27 | 27cm <b>\$124<sup>95</sup></b> |
| REF 3590-30 | 30cm <b>\$129</b> 95           |



#### Angled

| REF 3595-21 | 21cm | \$ <b>114</b> 95 |
|-------------|------|------------------|
| REF 3595-24 | 24cm | \$ <b>119</b> 95 |
| REF 3595-27 | 27cm | \$ <b>124</b> 95 |
| REF 3595-30 | 30cm | \$ <b>129</b> 95 |

# **HYSTERECTOMY SCISSORS**



# Titanium Rainbow<sup>™</sup> Scissors

#### **Titanium Coated Hysterectomy Scissors**



#### **Slight Curve**

| <b>REF TT680-20</b> | 20cm | \$ <b>97</b> <sup>95</sup> |
|---------------------|------|----------------------------|
| <b>REF TT680-23</b> | 23cm | \$ <b>99</b> <sup>95</sup> |



#### **Full Curve**

| <b>REF TT690-20</b> | 20cm | \$ <b>97</b> <sup>95</sup> |
|---------------------|------|----------------------------|
| REF TT690-23        | 23cm | \$ <b>99</b> <sup>95</sup> |



#### Curve

| <b>REF TT685-20</b> | 20cm | \$ <b>97</b> <sup>95</sup> |
|---------------------|------|----------------------------|
| <b>REF TT685-23</b> | 23cm | \$ <b>99</b> <sup>95</sup> |



#### Angled

| <b>REF TT695-20</b> | 20cm | \$ <b>97</b> <sup>95</sup> |
|---------------------|------|----------------------------|
| <b>REF</b> TT695-23 | 23cm | \$ <b>99</b> <sup>95</sup> |

#### Zephyr<sup>™</sup> Scissors

#### Stainless Steel Hysterectomy Scissors



#### **Curve Scissors**

| REF 3685-20 | 20cm | \$ <b>97</b> <sup>95</sup> |
|-------------|------|----------------------------|
| REF 3685-23 | 23cm | \$ <b>99</b> <sup>95</sup> |



#### **Angled Scissors**

| REF 3695-20 | 20cm | \$ <b>97</b> <sup>95</sup> |
|-------------|------|----------------------------|
| REF 3695-23 | 23cm | \$ <b>99</b> <sup>95</sup> |



#### **Slight Curve Scissors**

| REF 3680-20 | 20cm | \$ <b>97</b> <sup>95</sup> |
|-------------|------|----------------------------|
| REF 3680-23 | 23cm | \$ <b>99</b> <sup>95</sup> |



#### **Full Curve Scissors**

| REF 3690-20 | 20cm | \$ <b>97</b> <sup>95</sup> |
|-------------|------|----------------------------|
| REF 3690-23 | 23cm | \$ <b>99</b> <sup>95</sup> |

# **SCISSORS**

### Iris Scissors

#### Straight

REF \$3816 11cm Length......\$2895

#### Curved

REF \$3814 11cm Length......\$3495

# -----

#### Angled REF S3812

11cm Length......\$3995



# **Laminaria Insertion Forceps**

Specifically designed for insertion of laminaria. **REF** 1600 ......<sup>\$7695</sup>

# **NEEDLE HOLDERS**

### Webster Smooth Needle Holder

#### Smooth jaw

REF NHW3400

13cm Length......\$3295

# Mayo-Hegar Needle Holder

#### Serrated jaw

| REF NHMH3410-14 | 14cm Length <b>\$29</b> 95 |
|-----------------|----------------------------|
| REF NHMH3410-16 | 16cm Length <b>\$31</b> 95 |
| REF NHMH3410-18 | 18cm Length <b>\$35</b> 95 |
| REF NHMH3410-20 | 20cm Length\$3995          |

### **Crile-Wood Needle Holder**

#### Serrated jaw

| REF NHCW3420-15 | 15cm Length <b>\$31</b> 95 |
|-----------------|----------------------------|
| REF NHCW3420-18 | 18cm Length <b>\$35</b> 95 |
| REF NHCW3420-20 | 20cm Length <b>\$39</b> 95 |
| REF NHCW3420-23 | 23cm Length\$ <b>45</b> 95 |

Same Day Shipping for All In-Stock Orders Received by 2PM CDT

### **Littaur Stitch Scissors**

**REF S3810** 14

14cm Length......\$2995





#### **Electrode Extender**

#### Effectively Doubles the Shaft Length of Your Electrode

STERLE EO (Single-Use 🕲) Not made with natural rubber latex. **REF** 12-0000 13cm Length Box of 5 .....\$4295 ACCESS DEEP CAVITIES AND OTHER **HARD-TO-REACH SITES** 

# Electrode Handpiece, Hand Activated with VL Plug

#### For use with standard removable electrodes

 STERLE EO
 Single-Use (Single-Use (Si

**REF 755VL** Box of 12..... **\$105**<sup>95</sup>



# Electrode Handpiece, Foot Activated

#### For use with standard removable electrodes

STERLE EO Single-Use Single Sterle St


Not made with natural rubber lates

10 foot cable

**REF** 950

Box of 12......\$6895






# **H-Series Loop Electrodes**

Precision Depth and Width Control








# **GYNEX®** Cone Electrodes

STERLE EO Single-Use S Not made with natural rubber latex

#### The Cone Electrode is a LLETZ excisor for precise cone biopsy specimens


Cone biopsy electrodes allow for smooth rotation, facilitating the removal of precise, cone-shaped specimens. The tungsten linear wires are 0.20mm thick to help reduce backward drag. Molded shaft and stop arm are completely insulated to protect the endocervical canal and underlying vaginal structures.


Each electrode is packaged in a molded blister tray which keeps the tip true to shape and protects from breakage.



**Triangle Cone Electrode** 


### **Contoured Cone Electrode**






New

| REF 14-1020 | 10mm Radius x 20mm Height  |
|-------------|----------------------------|
| Box of 5    | \$ <b>59</b> <sup>95</sup> |





# **Electrosurgical Patient Dispersive Pads**

for LLETZ Procedures

NON-STERILE Single-Use Single States.








# Grounding Pad with Cord, Split (REM)





# Reach Kit<sup>M</sup> For Practical Paracervical Blocks

STERLE EO Single-Use (Single-Use (Single-U

The Gynex Reach Kit<sup>™</sup> allows you to **REACH** the injection site with ease. The kit contains everything you need to perform a paracervical block with the option of two convenient lengths. The Gynex Reach Kit<sup>™</sup> is the perfect solution for those looking to reduce cost without compromise. The Reach Kit<sup>™</sup> Cervical Block Procedure Kit is available in two different lengths.

#### 3" Needle Extender

Syringe Hub to Needle Tip: 10.1cm

 REF 5130
 Box of 6 Kits
 \$34<sup>95</sup>

#### 4" Needle Extender

Syringe Hub to Needle Tip:12.1cm

**REF 5140** Box of 6 Kits ...... **\$34**<sup>95</sup>

#### **Each Kit Contains:**


- Syringe, 10ml
- Draw-Up Needle, 18 GA
- Needle Extender, 3" or 4"
- Injection Needle, 0.5" 27 GA

### Extended Needle Adapters, Luer Hub, REUSABLE

Extends the reach of hypodermic needles

Stainless Steel, Available: 3-6 inch lengths





### Cartridge Syringe, Screw Hub, REUSABLE

NOTE: This does not work with the Extended Needle Adapter shown above.

REF 5210 \$52<sup>95</sup>

SALE

## Filter, Prefilter and Reducer

for Smoke Evacuation Systems

## Sterile Vacuum Tubing Sets, Single-Use

Fits many popular smoke evacuation systems that use 1/2" OD. The 9'+ tubing is connected to DSE tubing which connects to the speculum as shown in the illustration above.

| STERLE                             | ΕO | Single-Use (2) |
|------------------------------------|----|----------------|
| Not made with natural rubber latex |    |                |

**REF 8200** 12 Sets.....**\$198**<sup>95</sup>

### Prefilter Reducer, Single-Use

Fits into top of Prefilter (**REF** 8210), and adapts to .375" ID tubing (**REF** 8200 Tubing Set)

#### (NON-STERILE) Single-Use (8)

- Not made with natural rubber latex  ${\scriptstyle \blacksquare}$
- Bottom opening is 1.33" OD
- Top opening is .375" OD
- Fits CooperSurgical 6080 Smoke Evacuator, with 8210 Prefilter

REF 8230 Qty 10 ...... \$52<sup>95</sup>

### Prefilter, Single-Use

### NON-STERILE Single-Use (8)

#### Not made with natural rubber latex $\blacksquare$

- Fits CooperSurgical 6080 Smoke Evacuator
- Top opening is 1.33" ID
- Non-combustible

REF 8210

Qty 10 ..... \$189<sup>95</sup>

### ULPA Filter, Single-Use

#### NON-STERILE Single-Use (8)

Not made with natural rubber latex  $\blacksquare$ 

- Maximum Life 90 Days
- Fits CooperSurgical 6080 Smoke Evacuator



Perfect match for the competitor's smoke evacuator at a far more reasonable price!



CONNECTS WITH


CONNECTS WITH





## **Tubing and Connectors**

for LLETZ Procedures



## **STERILIZATION NOTICE** Please read before ordering Gynex LLETZ coated instruments

We find that most autoclaves are default set at **270**°. Gynex LLETZ coating must be autoclaved at **250**°. Please lower autoclave setting to this temperature to avoid damage to the coating. Bubbling can occur at a higher temperature. This also voids the warranty on these coated instruments.

#### Please use the following instructions:

- Use only distilled water in the autoclave.
- Assure that any towels the instruments are laid out on are free of residual soap, and are essentially pH neutral.
- Autoclave instruments constructed of like materials at the same time. E.g., autoclave only GYNEX LLETZ-Coated instruments together.
- The recommended sterilization cycle is **saturated steam** at **250°** for **25–30 minutes**. Use of this cycle will produce a Sterility Assurance Level (SAL) of at least 10<sup>-6</sup>.
- Do not allow instruments to touch each other during the autoclaving process, and assure that all instruments are **completely dry** when removed from the autoclave.

## Sterrad Sterilization Systems, can be used for Gynex LLETZ coated instruments.

Sterrad 100NX system has been established by Advanced Sterilization Products as the testing model for material functionality, representing all Sterrad System models. This can be applied to Sterrad NX, 100S, 50 and 200 systems and will be at least equivalent. Follow the manufacturer's instructions for operation and loading of sterilized items. Coating is compatible with at least 30 sterilization cycles with the 100NX System being the most challenging system on material functionality, although 60 cycles is the average compatibility. After 30 cycles, surface bubbling may be visible, but this will not impact the functionality of the product. If cracking of the coating is observed, do not use the instrument.

## Sponge Forceps, Single-Use

NON-STERILE Single-Use Not made with natural rubber latex Length: 24cm

 REF 3207
 12-Pack......\$39<sup>95</sup>

 REF 3214
 24-Pack ......\$73<sup>95</sup>



## **TOPICAL ANESTHETIC**

### HurriCaine<sup>®</sup> Topical Anesthetic Snap-n-Go<sup>™</sup> Swabs

Box of 72

- Fast onset 20 seconds
- Dose .15ml

See product details from manufacturer. Beutlich® Pharmaceuticals, LLC

REF 0569 \$53<sup>95</sup>

## HurriCaine<sup>®</sup> Topical Anesthetic Spray Kit

2oz Spray Container

- Provides 200 one second applications
- Includes 200 Single-Use extension tubes

See product details from manufacturer. Beutlich® Pharmaceuticals, LLC



REF 0183 \$72<sup>95</sup>





# LLETZ INSTRUMENTS

### **LLETZ-Coated Teneculum**

Extremely gentle tissue pickup • Coated Stainless Steel • Length: 25cm



5% Discount on Prepaid Orders by Check or ACH payment (offer excludes microscopes, colposcopes and/or extended warranties).


## **4-Way Expander Vaginal Speculum**

### for LLETZ Procedures

The 4-Way Expander Vaginal Speculum designed for LLETZ procedures with Gynex's autoclavable lavender coating. It provides a superior view with one instrument as compared to a standard speculum and lateral retractor. The built-in side blades are erected on demand to enhance the view of the cervix as needed.

#### Features






Blade Reference Chart: Shown Actual Size, Scale 1:1

## **LLETZ Open-Sided Speculum**

Left and Right Opening

Blade Reference Chart: Shown Actual Size, Scale 1:1



<sup>5%</sup> Discount on Prepaid Orders by Check or ACH payment (offer excludes microscopes, colposcopes and/or extended warranties).

## **Sullivan Vaginal Speculum**

for LLETZ Procedures


### **Features**

- Unique C-Arm maximizes your view
- Easy one handed open and release feature
- Locks firmly into place

#### ideal for bariatric patients

Blade Reference Chart: Shown Actual Size, Scale 1:1





Ideal to use with a Lateral Vaginal Retractor as there are no side obstructions like most specula









#### **Pederson & Graves** for LLETZ Procedures LLETZ Coated OK to Opening Shown Actual Size, Autoclave Scale 1:1 Read Notice on Page 74 Before Ordering **WideView** 40mm Blade Reference Chart: Shown Actual Size, Scale 1:1 Graves Pederson 180mm 175mm GYNEX 170mm 165mm 160mm 155mm 150mm 145mm 140mm 135mm 130mm 125mm 120mm 115mm 110mm 105mm 100mm bariatric bariatric 95mm 90mm 85mm 80mm 75mm 70mm 65mm 60mm **WideView** WideView WideView WideView WideView 55mm **REF 2845** REF 2847 50mm 45mm **FullView FullView FullView FullView FullView** 40mm REF 2883 **REF 2857 REF** 2887 35mm 30mm \$146<sup>95</sup> \$156<sup>95</sup> \$156<sup>95</sup> \$156<sup>95</sup> \$146<sup>95</sup> 25mm 20mm 15mm 10mm 5mm 80 0mm








## **Sullivan Stainless Steel Vaginal Speculum**


### **Features**


- Unique C-Arm maximizes your view
- Easy one handed open and release feature
- Locks firmly into place

### ideal for bariatric patients












C-Arm is designed to parallel the contour of the patient's leg

No side obstruction like most specula




## **Lateral Vaginal Retractor**

Control side wall collapse and get more cervical access with a GYNEX® Stainless Lateral Vaginal Retractor. For unparalleled ease of use and variety of practical sizes, nothing compares with GYNEX® Stainless Lateral Vaginal Retractor.

- Compatible with GYNEX<sup>®</sup> Wide-View<sup>™</sup> and FullView<sup>™</sup> Specula
- Lateral vaginal retractors fit inside the speculum providing greater visibility and access to the cervix





SALE

Blade Reference Chart: Shown Actual Size, Scale 1:1

## **Auvard Weighted Speculum**

Single-end vaginal speculum designed to be used during vaginal surgery (hysterectomy, dilation and curettage) with the patient in the lithotomy position. The weight holds the speculum in place to provide an open view.

Fixed Weight: 1.4 kg Size:105mm long x 38mm wide

REF SA2200

\$**299**<sup>95</sup>

## **Guttman Speculum** with Special Blades

This indispensable speculum operates without upper blades, permitting total visualization of the cervix & vagina. Two separators aid in the opening and stabilizing the vagina.

### Medium


| REF 2310 | Top Separator: 82mm |
|----------|---------------------|
| Large    |                     |
| REF 2312 | Top Separator: 92mm |




00.00

## Greenberg Hysteroscopy Speculum

The Greenberg Hysteroscopy Speculum is a modified traditional open-sided Graves speculum with shorter blades and an open side allowing removal of instruments while they are still in the vagina. This speculum brings the patient's cervix closer to the provider and reduces the length of unsupported hysteroscope shaft entering the cervical canal.





James Greenberg, MD, Boston, MA

| 90mm |           |  |
|------|-----------|--|
| 85mm |           |  |
| 80mm | Greenberg |  |
| 75mm |           |  |
| 70mm |           |  |
| 65mm |           |  |
| 60mm |           |  |
| 55mm |           |  |
| 50mm |           |  |
| 45mm |           |  |
| 40mm |           |  |
| 35mm |           |  |
| 30mm |           |  |
| 25mm |           |  |
| 20mm |           |  |
| 15mm |           |  |
| 10mm |           |  |
| 5mm  |           |  |
| 0mm  |           |  |
|      |           |  |


## **Open-Sided Speculum, STAINLESS**



**Note**: Blade length is measured from top center of back of blade to end of distal tip. In most cases, the bottom speculum blade will protrude and appear to be longer. Mechanical assembly of a speculum simply places the bottom blade in a more forward position. Although the two blades are normally similar in length, they may vary a few millimeters from lot to lot.

## Titanium Rainbow<sup>™</sup> Speculum

Exclusive Titanium Rainbow<sup>™</sup> coating extends the life of instruments over traditional stainless-steel. The titanium coated speculum offers increased resistance to corrosion, resilience against impact and low friction coefficient. Titanium Rainbow<sup>™</sup> is available in two viewing widths for Pederson and Graves.



## **GYNEX® ULTRASONIC FETAL DOPPLER**

Add this excellent Fetal Doppler to your practice!

## Fetal Doppler Model D2500

- Highly sensitive ultrasound
- Detects heart rates of fetus as early as the 12th week
- Automatic shutdown after 1 minute without operation
- Built-in loudspeaker
- Rechargeable battery, 1180 mah lithium
- Colorful display
- 2.5MHz Ultrasound Frequency
- FHR value, FHR waveform volume
- Includes battery charger and earphones

Dimension: 215mm x 200mm x 66mm

REF D2500 ......\$295<sup>95</sup>



## **AVAILABLE ACCESSORIES**

### **Thermasonic Gel Warmer**

- Maintains constant temperature
- Low electrical demand
- Heat indicator lamp
- Adjustable thermostat

REF D82-01 ......\$290<sup>95</sup>



### **Doppler Gel**

- Acoustically correct for the broad range of frequencies used
- Completely aqueous, will not stain clothing or damage transducers
- 250ml bottle

REF D-GEL Box of 6 ......\$39<sup>95</sup>



## pHem-Alert<sup>®</sup> Vaginal pH Test

(NON-STERILE) (Single-Use (2)) Not made with natural rubber latex =


pHem-Alert<sup>®</sup> test measures vaginal pH, and is intended for use by women who have any of the following vaginal symptoms: Itching, Burning, Unpleasant Odor, Unusual Discharge

### **KEY ADVANTAGES**

- Each test securely packed in a procedure blister tray
- Convenient, unique, patented applicator "key" with attached pH paper
- Easy to read pH chart included with every test set

**REF pH-99** Box of 12 ......<sup>\$</sup>**67**<sup>95</sup>

## The Original Vaginal pH Test Designed for Accurate & Comprehensive Results!



#### pH 4.5

- Normal pH
- Predominant lactobacilli
- Possible yeast overgrowth
- Normal serum estradiol

#### pH 5.0

- Low grade bacterial vaginitis
- Borderline serum estradiol

#### pH 5.5

- Bacterial vaginitis
- Low serum estradiol

#### pH 6.0

- Bacterial vaginitis
- Low serum estradiol

#### pH 6.5

- Bacterial vaginitis
- Very low serum estradiol

#### pH 7.0

Low to no serum estradiol

#### pH 7.5

Low to no serum estradiol

## **INSTRUMENTS CARE AND CLEANING**

## Touch-Guard® Hot Beads

### **Sterilization Beads**

- Touch-Guard<sup>®</sup> Hot Beads provide a safe, effective and convenient way to sterilize surgical instruments.
- No need for potentially dangerous gases, flames, or chemicals to achieve sterilization.
- The specially cast "no-lead" glass beads transmit heat properly and remain intact at high temperatures.
- Simply plunge metal instruments into the hot glass beads for instant sterilization.
- Touch-Guard<sup>®</sup> Hot Beads are free of pathogens and microbial contaminants.
- Pre-measured containers make it easy to use the right amount of beads for each sterilization cycle.



## Touch-Guard<sup>®</sup> Stain-Ex<sup>™</sup>

### **Biological & Chemical Stain Remover**

- Touch-Guard<sup>®</sup> Stain-Ex<sup>™</sup> is a fast and effective biological and chemical stain remover.
- It removes a variety of stains, including lodinebased stains, Lugol's stain, and medical/ histological stains.
- It can be used on a range of surfaces, including hands, glassware, stainless steel, plastic, counters, floors, lab coats, laboratory equipment, and hospital linen.
- Touch-Guard<sup>®</sup> Stain-Ex<sup>™</sup> works great on LLETZcoated instruments.
- Application is easy, and stains are removed quickly (in a matter of seconds).





## **Touch-Guard® Low Suds**

### Liquid Instrument Detergent Concentrate

Touch-Guard<sup>®</sup> Low Suds Liquid Instrument Detergent Concentrate is used for effectively and safely cleaning all surgical instruments, plastic, glassware, and rubber tubing. Formulated for use in automatic washers, ultrasonics, washersterilizers, and manual cleaning.

- Neutral pH detergent with powerful sequestering and emulsifying agents
- Fast wetting, odorless, non-sudsing, maintains high detergency
- Removes protein, blood, tissue, fat and any organic soils
- No residue no pyrogens; leaves a surgically clean surface
- Exceptional cleaning in hard, soft, cold or hot water
- Biodegradable; contains no phosphates or NTA
- Free rinsing, non-enzymatic formulation
- Safe for ophthalmic cleaning

REF TG-200







## **Touch-Guard® Lubricant**

### Instrument Lubricant and Rust Inhibitor

Touch-Guard<sup>®</sup> Instrument Lubricant and Rust Inhibitor protects and extends the life of all your surgical instruments. When used daily as a rinse after cleaning, it will keep your instruments moving parts operating smoothly and remain effective even after autoclaving, without any oily residue.

- Prevents rust, spots, stains and corrosion
- · Restores articulation to sticky box-locks and joints
- Withstands high temperatures ad autoclaving
- Steam penetrable, sterizable and ETO sterizable
- Non-oily, non sticky, non-toxic and silicone free



**REF TG-100** 16.9 fl. oz. (500 ml) ......**\$16**<sup>95</sup>

FREE

SAMPLE

To redeem your free sample, contact us (888) 486-4644 Toruch-Guarden Toruch-Guarden Instrument Lubricati and Russ Inhibitor We tai Mission Source Source Source Toruch Control Source Source Source Mission Source Source Source Source Source Mission Source Source Source Source Source Source Mission Source So

man and represents the reducting sping wat the memory of and summerics that its cosing constitutions

> To Rearder: Prove 1 and an alle Fax 1.403,001,011 and environmental

GYNEX

## **INSTRUMENT REPAIR AND SHARPENING**

Gynex<sup>®</sup> allows you to extend the life of your instruments, deliver optimal performance and better manage your surgical asset investment. All sharpening and repairs also include rust/stain removal, cleaning, minor polishing and lubrication.

Gynex<sup>®</sup> sharpens and repairs most medical instruments. If the instrument that you would like to have sharpened or repaired is not listed here, or you have any questions, please contact us for more information.

Please note that due to metal, age and condition, some instruments may not be repairable or sharpened.

### Pricing is as follows:

Gynex® Biopsy Punch ......<sup>\$</sup>65<sup>00</sup> Other Brand Biopsy Punch ......<sup>\$</sup>75<sup>00</sup> Gynex® Curette ......<sup>\$</sup>40<sup>00</sup> Other Brand Curette ......<sup>\$</sup>50<sup>00</sup>

**Toll Free: (888) 486-4644** (8:00AM-5:30PM CDT)

#### **Selling Notice:**

All products are limited to purchase by licensed medical practitioners, physician and medical facilities. In most cases, we do not sell to the general public.

#### Warranty:

All hand-held instruments in this catalog are warranted to be free of defects in manufacturing materials and workmanship for two years from date of purchase. Customers who feel an instrument is defective must return it to GYNEX for evaluation. Products found to be defective will either be repaired or replaced with the same or comparable product, at the discretion of GYNEX Customer Service. Warranty on any hand-held GYNEX instrument is voided if instrument is repaired or sharpened by anyone other than GY-NEX. All GYNEX biopsy punches in this catalog will be sharpened free of charge for one year from date of purchase (if enough metal is remaining on tip to allow sharpening). Warranties do not cover general maintenance, cleaning or sharpening needs that result from normal use, except as noted above regarding biopsy punch sharpening. Damages caused by misuse, improper care, cleaning, chemicals, or loss and theft are not covered by this warranty. Shipping charges are not covered for free sharpening warranty nor extended equipment warranty.

#### **Return Policy:**

Customer must obtain а Returned Merchandise Authorization (RMA) number from GYNEX customer service at (888) 486-4644. An itemized list of returned products must accompany all returns. Include the Returned Merchandise Authorization (RMA) number and explanation of why the product is being returned. Returned goods MUST be packaged in a shipping container equivalent to that in which it was received. Returned goods MUST be shipped at customer's expense by Parcel Post, Airborne, FedEx or UPS, and must be insured for the purchase price of the contents. Returned goods received by GYNEX in any form less than that described above will be refused by GYNEX and return shipped to the customer. The following items may be returned:

- Product with RMA number
- All returns must be made with in 30 days of the item receipt (except for warrantv).
- Product with unbroken seals, with labels intact, in original container and packaging
- Unused instruments which have not been cleaned or sterilized
- Single-Use<sup>®</sup> product (sterile) in original, unopened package
- Product damaged in shipment\*
- Product that has at least 12 months remaining shelf life as determined by the expiration date
- Instruments damaged, broken or not working properly due to defects in material or workmanship
- \*Note: Product returned due to error or damage in shipping must be reported to customer service within 72 hours of receipt. Damaged product must be accompanied by a written explanation noting damage.

## The following items may not be returned:

- Product without RMA number
- Product that has been improperly handled or stored
- Product that is marked, coded, dated, damaged, soiled or adulterated in any way (except standard cleaning and sterilization procedures for reusable instruments)
- Product damaged by or exposed to heat, fire, smoke, high humidity or water (except standard cleaning and sterilization procedures for reusable instruments)

#### **Restocking Charge:**

A 15% restocking charge applies to all returned product, unless the return is due to error in shipment, product damaged in shipment, or product returned for repair or replacement due to defect in material or workmanship.

#### **Disclaimer:**

All prices in this catalog are subject to change without notice. All sizes indicated may vary slightly, due to differences between manufacturers. In order to carry sufficient inventory, it is necessary for GY-NEX to purchase from different manufacturers. As a result, final products may not be the exact sizes as shown in charts, drawings or photos. The advertisements and information in this catalog may contain typographical errors or inaccuracies and may not be complete or current. GYNEX therefore reserves the right to correct any errors, inaccuracies, or omissions and change or update information at any time without prior notice. Please note that such errors, inaccuracies, or omissions may relate to pricing or availability. We reserve the right to refuse to honor offers or terms, which contain error or inaccuracies.

#### **Trademark Notices:**

Apex, Ascona, Ceramix, Century Uterine Manipulator/Injector. Clear-Flue, C-Safe, Direct-Access, EasyView, HexLoc, i-Series Titanium LiteTouch, LongView, Rainbow, LVR, NoSlip, FullView, New GYNeration, Opal Rainbow, Reach, Reach Extended Needle, RotoTip, Stain-Ex, Titanium Rainbow, Touch-Guard, Villa, WideView, Willows and Zephyr are trademarks of GYNEX Corporation. Dennis Dilators and Galileo are registered trademarks of GYNEX Corporation. . Use of these product names is limited to those with written authorization by an officer of **GYNEX Corporation of North Richland** Hills, Texas. ConMed and Sabre are registered trademarks of ConMed, Inc. UMI is a trademark of Thomas Medical, Inc. CooperSurgical, Shur-Lock and Euro-Med are registered trademarks of Cooper Companies. HurriCaine is a registered trademark and Snap -n- Go is a trademark of Beutlich LP, Pharmaceuticals. pHem-Alert is registered trademark of FemTek, LLC. Suresample is a registered trademark, and Portex is a trademark of Smiths Group. Z-Sampler is a trademark of BEI Medical Systems Company Inc. BEI is a registered trademark of BEI Medical Systems Company Inc. Pipelle is a registered trademark of Unimar, Inc. Utah Medical Products and Finesse are registered trademarks of Utah

Medical Products, Inc. Endocell is a trademark of Wallach Surgical Devices, Inc., and Wallach is a registered trademark of Wallach Surgical Devices, Inc. Teflon is a registered trademark of DuPont.

#### Force Majeure:

Neither party shall be liable for any failure or delay in delivery of goods, or in otherwise fulfilling its obligations with respect to the sale of goods to the other which results from any cause beyond its reasonable control, including disruption or delays in a common carrier's transportation schedule, compliance with any regulations, orders, acts, instructions, or priority requests of any national, state or municipal government, or any department or agency thereof, or any civil or military authority, or from acts of God, fires, or other casualty of accidents, strikes, lockouts, factory shutdowns, embargoes, riots, or other disorders, delays or shortages in transportation, or inability to obtain in sufficient quantity fuel, power, labor, manufacturing facilities or materials or other supplies from the usual sources of GYNEX. Not with standing the foregoing, in the event a party is unable to perform as a result of such force majeure event for a period exceeding sixty (60) days, the other party may, at its option, terminate its purchase order or commitment to sell the goods and neither party shall have any further liability hereunder other than for past services rendered and goods sold.

#### **Official Notice**

ALL products are shipped from GYNEX state-of-the-art warehouse in North Richland Hills, Texas


## **INDEX**

| ACETIC ACID2  |
|---|
| Alligator Grapser, Gynex <sup>®</sup>                             |
|   |
| Ascona™ IUD String Snare  |
| Accu-Sound <sup>®</sup>   |
| BARIATRIC INSTRUMENTS   |
| Biopsy Punches6-14  |
| Tissue Forceps  |
| Specula, LLETZ Coated80-81  |
| Specula, Stainless Steel82-83                                     |
| BIOPSY PUNCHES  |
| Amethyst  |
| Apex <sup>™</sup>   |
| Ceramix <sup>™</sup> 6  |
| Gold  |
| Kevorkian Ring Handle 14  |
| Medium-Tischler 11  |
| Mini-Tischler 10  |
| Opal Rainbow™7  |
| Rose Gold   |
| Sapphire  |
| ROTOTip <sup>™</sup> Biopsy Punches12-13                          |
| Sullivan Ring Handle 14   |
| CARE AND CLEANING   |
| Hot Glass Beads Sterilizer - Touch-Guard®92                       |
| Instrument Lubrication - Touch-Guard <sup>®</sup>                 |
| Instrument Repair and Sharpening                                  |
| Low Suds Detergent - Touch-Guard <sup>®</sup>                     |
| Stain-Ex <sup>™</sup> Stain Remover - Touch-Guard <sup>®</sup> 92 |
| COLPOSCOPES   |
| Galileo®46-51   |
| Accessories for Galileo <sup>®</sup>                              |
| CO Series   |
| DILATORS  |
| All Uterine Dilators  |
| Dennis Dilators, Single-Use (8)                                   |
| Dennis Dilators, Reusable   |
| Hank Dilators, Uterine  |
| Hegar Dilators, Uterine   |
| Mini Dilators, Uterine  |
| Mini Dilators Tips, 4mm and 5mm                                   |
| Pratt Dilators, Uterine   |
| Small Dilators, Uterine   |
| Soft Dilators, Single-Use (8)                                     |
| Soft Teflon Dilators, Reusable                                    |
| Doppler Gel90 ELECTRODES  |
| Ball Electrode  |
| Contoured Loop Electrodes   |
| Electrode Extender  |
|   |
| Electrode Handpieces65<br>H-Series Loop Electrodes67              |
| Needle Electrode  |
| Radius Loop Electrode   |
| Triangle Loop Electrodes  |
| 11 augre Loop Liebu dues  |
| Square Loop Electrodes 67   |
| Square Loop Electrodes  |

| ENDOCERVICAL   |
|--|
| Curettes - Rainbow™and Stainless   |
| Kevorkian Curette, Single-Use  |
| Lickerish Spiral Curettes 17   |
| Specula 19   |
| ENDOMETRIAL CURETTE AND SAMPLER  |
| SURESAMPLE® Select1 Curette42-43   |
| SURESAMPLE® Sampler44-45   |
| FORCEPS  |
| Foerster Forceps, Stainless Steel  |
| Laminaria Insertion Forceps64  |
| Sponge Forceps, Single-Use74   |
| Ovum Forceps 20-21   |
| FETAL DOPPLER - GYNEX®90   |
| GRASPER  |
| Alligator Grapser - GYNEX®28   |
| Serpentine Grasper28   |
| HOOKS  |
| LLETZ Coated75   |
| Stainless Steel  |
| IUD INSTRUMENTS  |
| Ascona <sup>™</sup> IUD String Snare                                       |
| Ceramix <sup>™</sup> Coated String Scissors                                |
| Crochet Type IUD Removers  |
| GYNEX® Alligator Grasper   |
| IUD Insertion Kit, Single-Use -GYNEX <sup>®</sup>                          |
| IUD String Scissors  |
| Serpentine Grasper IUD Remover   |
| Single-Use IUD String Scissors   |
| Villa <sup>™</sup> Hormonal Dissecting Forceps                             |
| Villa™ Hormonal Implant Remover27<br>Willows™ IUD String Snare, Single-Use |
| LATERAL VAGINAL RETRACTORS   |
| LATERAL VAUIVAL RETRACTORS   |
| Stainless Steel  |
| LLETZ COATED INSTRUMENTS   |
| All Items  |
| Air reens  |
| LVR  |
| Speculum   |
| Sponge Forceps, Single-Use ®74   |
| MEDICAL EXAM LIGHT - Galileo®  |
| Medical Exam Light - Accessories   |
| Medical Exam Light Protective Sheath                                       |
| Microscope - Galileo®53  |
| MONSEL'S PASTE   |
| NEEDLE ADAPTER   |
| NEEDLE HOLDERS   |
| Webster64  |
| Mayo-Hegar64   |
| Crile-Wood64   |
| OVUM FORCEPS   |
| Bierer   |
| Blumenthal21   |
| Jacobson21   |
| Saenger  |
| Sopher Original  |

| PATEL <sup>™</sup> MULTI-USE FORCEPS                        | ļ |
|---|---|
| PUNCHES   |   |
| Keyes   |   |
| Keyes Accessories22   |   |
| REACH KIT <sup>™</sup> Cervical Block Procedure Tray2,70    |   |
| REPAIRS & SHARPENING SERVICE                                | 3 |
| SCISSORS  |   |
| IUD String25  |   |
| Titanium Rainbow™63   |   |
| Zephyr™Stainless Steel63                                    |   |
| Iris Stainless Steel64                                      |   |
| Littaur Stitch Scissors64                                   | ļ |
| SMOKE EVACUATION  |   |
| Prefilter, Single-Use <sup>®</sup> 72                       |   |
| Prefilter Reducer72   | 2 |
| ULPA Filter72   |   |
| Tubing73  | 3 |
| SPECULA   |   |
| LLETZ Coated  |   |
| 4-Way Expander76  |   |
| Graves80-81   |   |
| Open-Sided (left/right opening)77                           |   |
| Pedersons80   |   |
| Sullivan78-79   | ) |
| Stainless Steel   |   |
| Auvard86  | j |
| Graves82-83   | 3 |
| Greenberg87   | / |
| Guttman87   |   |
| Open-Sided (left/right opening)88                           | 3 |
| Pederson82  | 2 |
| Sullivan84-85   | j |
| Titanium Rainbow™   |   |
| Graves89  | J |
| Pederson89  |   |
| STERILIZATION NOTICE74                                      |   |
| Syringe, Cartridge®71                                       |   |
| TABLE OF CONTENTS 1   |   |
| TENACULA  |   |
| Braun   | 3 |
| Schroeder18   | 3 |
| Skene19   | J |
| Teale19   | J |
| TOPICAL ANESTHETIC - HurriCaine®                            |   |
| Spray Kit74   | ł |
| Swabs74   |   |
| TRANSIT MAP96   |   |
| UTERINE MANIPULATOR / INJECTOR - CenturyUMI <sup>™</sup> 41 |   |
| UTERINE SOUNDS  |   |
| Accu-Sound®31   |   |
| Sim Sound32   |   |
| Sim Sound • flex32  |   |
| VAGINAL PH TESTING - pHEM-ALERT®91                          |   |
| WORD BARTHOLIN CATHETER KIT24                               | ł |

## Transit Map



Fax: (425) 895-0115 E-mail: service@gynex.com Web: www.gynex.com

**Terms:** Net 30 Days on open approved accounts. 1.0% finance charge per month shall apply on all past due accounts. International accounts require prepayment. New accounts may require prepayment with first order.

**Shipping:** FOB North Richland Hills, Texas. Normal shipping via UPS or FedEX Ground, and charges for shipping and handling are added to the order total. Expedited delivery (Next Day, 2nd Day, 3rd Day) available for additional charge.

**Note:** Pre-paid shipping charges of \$17.50, do not apply to Alaska, Hawaii, Puerto Rico, the Virgin Islands and other U.S. Islands, international orders or orders including colposcopes and medical exam lights. An additional custom fee may apply for the Virgin Islands and other U.S. Islands, and international orders. A handling fee will apply to orders that use customer's own shipping account.

**Partial Shipments:** We reserve the right to make partial shipments against orders and the amounts for such partial shipments will become due and payable upon receipt unless otherwise agreed to in writing by GYNEX.

**5% Discounts:** Customers may deduct an additional 5% off their total order when payment is made by check or ACH at the time of purchase. Credit card payments and COD orders do not qualify for this discount. This discount does not apply to Microscopes, Colposcopes and/or extended warranties and some specials where products are already discounted. This discount cannot be combined with other discounts.

Warranties Returns and Disclaimers: All information concerning these GYNEX policies are detailed on page 94 of this catalog.

Missing Merchandise Policy: GYNEX personnel have taken great care to pack your order. Each LOT number indicated in the first column on the Packing List corresponds with and was scanned from the product packed in this box. Quantities were verified in this manner. Back orders are clearly marked, and do not have a LOT number indicated. If anything is missing from your package, the customer is responsible for notifying GYNEX within 24 hours of the customer's receipt of package.

PayPal

WE ACCEPT

VISA

| GYNE   | EX®                    | MAIL TO: GYNEX<br>P.O. BOX 822369<br>N. RICHLAND HILLS, TX 76182 |
|--|------------------------|--|
| MAKE A COPY OF THIS ORDER FORM, FILL IT OUT AN | ID SEND BY FAX OR MAIL | FAX TO: (425) 895-0115   |
| BILLING ADDRESS                                | <b>SHIP TO</b>         |  |
| Name   | Name                   |  |
| Address  | Address                |  |
| City State Zip                                 | CitySta                | ite Zip  |
| Email  | Fax                    |  |
| Phone ( )                                      | Contact Person         |  |

| REF | Qty. | DESCRIPTION | Price Each | TOTAL |
|-----|------|-------------|------------|-------|
|     |      |             |            |       |
|     |      |             |            |       |
|     |      |             |            |       |
|     |      |             |            |       |
|     |      |             |            |       |
|     |      |             |            |       |
|     |      |             |            |       |
|     |      |             |            |       |
|     |      |             |            |       |
|     |      |             |            |       |
|     |      |             |            |       |
|     |      |             |            |       |
|     |      |             |            |       |
|     |      |             |            |       |
|     |      |             |            |       |
|     |      |             |            | )     |

|                              | Total Amount Ordered   |         |
|------------------------------|--|---------|
| MAKE CHECKS PAYABLE TO:      | 5% Discount for Pre-Paid Orders ONLY   |         |
| GYNEX CORP.<br>PO BOX 822369 | Shipping and Handling for Pre-Paid Orders ONLY<br>GROUND SHIPPING TO CONTINENTAL US ONLY, EXCLUDING ALASKA<br>AND HAWAII. (Shipping charges apply for orders to Alaska and Hawaii) | \$17.50 |
| N. RICHLAND HILLS, TX 76182  | Sales Tax  |         |
|                              | Total  |         |



6631 Browning Drive, North Richland Hills, TX 76180

Tel: (888) 486-4644 • Fax: (425) 895-0115 service@gynex.com • www.gynex.com





SCAN to SHOP





(888) 486-4644 www.gynex.com

